

# INSTRUCTION MANUAL

## **Quick-DNA/RNA™ Miniprep Kit**

Catalog No. **D7001**

### **Highlights**

- Quick, 15 minute isolation and separation of DNA and RNA (~25 µg) from a wide range of sources using Zymo-Spin™ column technology.
- High-quality DNA/RNA eluted in ≥25 µl is ready for reverse transcription, microarray, sequencing, etc.

### **Contents**

Product Contents .....	1
Product Specifications.....	1
Product Description.....	2
Reagent Preparation.....	3
Protocol.....	3, 4
Appendix.....	5
Ordering Information .....	6

Satisfaction of all Zymo Research products is guaranteed. If you are dissatisfied with this product please call 1-888-882-9682.

\* For assistance, contact us at [tech@zymoresearch.com](mailto:tech@zymoresearch.com).

#### Notes:

<sup>2</sup> For purification of DNA and RNA from whole blood, see the **Quick-DNA/RNA™ Plus Kits** (D7003, D7005).

This product is for research use only and should only be used by trained professionals. It is not for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

™ Trademarks of Zymo Research Corporation. RNAlater™ is a trademark of Ambion, Inc.

## Product Contents

<b>Quick-DNA/RNA™ Miniprep Kit (Kit Size)</b>	<b>D7001 (50 preps.)</b>	<b>Storage Temperature</b>
<b>DNA/RNA Lysis Buffer</b>	50 ml	Room Temp.
<b>DNA/RNA Prep Buffer</b>	50 ml	Room Temp.
<b>DNA/RNA Wash Buffer<sup>1</sup> (concentrate)</b>	24 ml	Room Temp.
<b>DNase/RNase-Free Water</b>	10 ml	Room Temp.
<b>Zymo-Spin™ IIC Columns</b>	50	Room Temp.
<b>Zymo-Spin™ IICG Columns</b>	50	Room Temp.
<b>Collection Tubes</b>	3x 50	Room Temp.
<b>Instruction Manual</b>	1	-

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

<sup>1</sup> Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **DNA/RNA Wash Buffer** concentrate before use.

## Specifications


- **Sample Types** – Cells, small amounts of *easy-to-lyse* tissue, buffy coat, buccal cells, plasma, serum, and other biological liquids. *Not compatible with whole blood.*<sup>2</sup>
- **Sample Size** – 10<sup>2</sup> to 5x10<sup>6</sup> cells in suspension or as tissue.
- **Yield** – DNA and RNA can be eluted into small volumes (≥25 µl) allowing for a highly concentrated sample. Maximum DNA/RNA binding capacity of the provided columns is ~25 µg.
- **Size Limits** – Capable of recovering genomic DNA up to and above 40 kb. In most instances, mitochondrial DNA and viral DNA (if present) will also be recovered. Total RNA ≥17 nucleotides can be recovered.
- **Purity** – High quality genomic DNA and total RNA ( $A_{260}/A_{280} > 1.8$ ,  $A_{260}/A_{230} > 1.8$ ) is recovered. Traces of DNA may be present in the eluted RNA fraction. Trace DNA can be removed by DNase digestion (page 5).
- **RNA Storage** – RNA is eluted with RNase-free water and can be stored at ≤-70 °C. The addition of RNase inhibitors is highly recommended for prolonged storage.
- **Required Equipment** – Microcentrifuge

### ZYMO RESEARCH CORP.


## Product Description

The **Quick-DNA/RNA™ Miniprep Kit** provides a quick method for the isolation of high quality genomic DNA and total RNA from small amounts of cells and tissue. The kit isolates *both* genomic DNA and a broad range of RNA species without the use of phenol. Small RNAs (*e.g.*, tRNAs, microRNAs) can be recovered following a simple adjustment within the RNA isolation protocol – *no extra steps are required!* Both DNA and RNA from up to  $5 \times 10^6$  cells can be eluted into volumes as little as 25  $\mu$ l in less than 15 minutes.

For **Assistance**, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail [tech@zymoresearch.com](mailto:tech@zymoresearch.com).



Genomic DNA (lane 1, 2) and total RNA (lane 3, 4) isolated from human epithelial cells (HCT 116) with the **Quick-DNA/RNA™ Miniprep Kit**.



PCR amplification of  $\beta$ -actin transcript (353 bp fragment shown) following DNA and RNA isolation from human epithelial cells (HCT 116) with the **Quick-DNA/RNA™ Miniprep Kit**: PCR positive control (DNA template; lane 1, 2, 3), PCR negative control (RNA template; lane 4, 5, 6), RT-PCR (lane 7, 8, 9).

### ZYMO RESEARCH CORP.

Ensure the RNA isolation procedure is performed in an RNase-free environment.

**Notes:**

<sup>1</sup> In order to lyse samples completely, the amount of the **DNA/RNA Lysis Buffer** can be adjusted (*i.e.*, more buffer can be added).

<sup>2</sup> The capacity of the **Zymo-Spin™ Column** is 800 µl. Columns can be reloaded to process volumes >800 µl.

<sup>3</sup> The maximum binding capacity of the **Zymo-Spin™ IIICG and IIC Column** is ~25 µg of DNA/RNA.

## **Reagent Preparation**

- ✓ Before starting, add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **DNA/RNA Wash Buffer** concentrate.

## **Protocol**

All centrifugation steps should be performed at 10,000 - 16,000 x *g* for 30 seconds unless specified. All steps should be performed at room temperature (20-30°C) unless specified.

### **1. Sample Preparation**

- A. Adherent Cells:** Cells can be lysed directly in the culture container by removing liquid medium and adding **DNA/RNA Lysis Buffer**<sup>1</sup> directly to the monolayer (*e.g.*, 400 µl for 10<sup>2</sup> to 5x10<sup>6</sup> cells). Remove cells from the culture surface by pipetting, scraping, *etc.* Proceed to *Step 2*.
  - B. Cells in Suspension:** Pellet the cells by gentle centrifugation (*e.g.*, 5 minutes at 500 x *g*). Remove the supernatant completely and resuspend the cell pellet in 400 µl **DNA/RNA Lysis Buffer**<sup>1</sup>. Vortex briefly. Proceed to *Step 2*.
  - C. Solid Tissue Samples:** Add 400 µl **DNA/RNA Lysis Buffer**<sup>1</sup> to fresh or frozen tissue (up to ~25 mg) and homogenize the sample (*e.g.*, using a Dounce or similar homogenizer). Proceed to *Step 2*.
  - D. Liquid Samples:** Add 3 volumes of **DNA/RNA Lysis Buffer**<sup>1</sup> for every volume of sample (*e.g.*, 300 µl buffer to 100 µl sample). Proceed to *Step 2*.
2. Transfer the sample from *Step 1* into a **Zymo-Spin™ IIICG Column**<sup>2,3</sup> in the **Collection Tube** and centrifuge at ≥12,000 x *g* for 1 minute.

*Save the flow-through for RNA purification and the column for DNA purification!*

**DNA Purification**

3. Transfer the **Zymo-Spin™ IICG Column** into a new **Collection Tube**.

4. Add 400 µl **DNA/RNA Prep Buffer** to the column and centrifuge at  $\geq 12,000 \times g$  for 1 minute. Discard the flow-through.
5. Add 700 µl **DNA/RNA Wash Buffer** to the column and centrifuge at  $\geq 12,000 \times g$  for 30 seconds. Discard the flow-through.
6. Add 400 µl **DNA/RNA Wash Buffer** and centrifuge the column for 2 minutes to ensure complete removal of the wash buffer. Carefully transfer the column into a clean microcentrifuge tube.

7. Add  $\geq 50$  µl **DNase/RNase-Free Water** directly to the column matrix and let stand 2 to 5 minutes at room temperature, then centrifuge at top speed for 30 seconds. The eluted DNA can be used immediately or stored at  $\leq -20^\circ\text{C}$ .

**RNA Purification**

3. Add an equal volume<sup>1</sup> of ethanol (95-100%) to the flow-through in the **Collection Tube**<sup>2</sup> from *Step 2* and mix well by pipetting. Then transfer the sample into a **Zymo-Spin™ IIC Column** in a **Collection Tube** and centrifuge at  $\geq 12,000 \times g$  for 1 minute. Discard the flow-through.<sup>3</sup>

7. Add  $\geq 25$  µl **DNase/RNase-Free Water** directly to the column matrix and let stand 1 minute at room temperature. Centrifuge at  $10,000 \times g$  for 30 seconds. The eluted RNA can be used immediately or stored at  $\leq -70^\circ\text{C}$ .

**Notes:**

<sup>1</sup> Alternatively, to isolate RNAs  $\geq 200$  nt, add  $\frac{1}{2}$  volume ethanol (95-100%) to the sample flow-through.

<sup>2</sup> Capacity of the **Collection Tube** is 2 ml.

<sup>3</sup> At this point, RNA samples can be in-column DNase I treated (page 5).

**Notes:**

<sup>1</sup> Prior to use, reconstitute the lyophilized **DNase I** as indicated on the vial. Store frozen aliquots.

*\* Unit definition - one unit increases the absorbance of a high molecular weight DNA solution at a rate of 0.001 A<sub>260</sub> units/min/ml of reaction mixture at 25°C.*

**Appendix A: In-Column DNase I Digestion**

The DNase I digestion procedure can be performed using **DNase I Set** (E1010).<sup>1</sup>  
All centrifugation steps should be performed at 10,000 –16,000 x g for 30 seconds unless specified.

1. Wash the column with 400 µl **DNA/RNA Wash Buffer** and centrifuge. Discard the flow-through.
2. Add 80 µl **DNase I Reaction Mix** (below) directly to the column matrix.

<b>DNase I</b>	5 µl (1 U/µl)*
<b>DNA Digestion Buffer</b>	75 µl

3. Incubate the column at room temperature (20-30°C) for 15 minutes.  
Continue with RNA Purification: Page 4, Step 4.

**Ordering Information**

<b>Product Description</b>	<b>Kit Size</b>	<b>Catalog No.</b>
<b>Quick-DNA/RNA™ Microprep Plus Kit</b>	50 Preps.	D7005
<b>Quick-DNA/RNA™ Miniprep Kit</b>	50 Preps.	D7001
<b>Quick-DNA/RNA™ Miniprep Plus Kit</b>	50 Preps.	D7003

<b>For Individual Sale</b>	<b>Amount</b>	<b>Catalog No.</b>
<b>DNA/RNA Lysis Buffer</b>	50 ml	D7001-1-50
<b>DNA/RNA Prep Buffer</b>	10 ml	D7010-2-10
	25 ml	D7010-2-25
	50 ml	D7010-2-50
<b>DNA/RNA Wash Buffer (concentrate)</b>	6 ml	D7010-3-6
	12 ml	D7010-3-12
	24 ml	D7010-3-24
<b>DNase/RNase-Free Water</b>	1 ml	W1001-1
	4 ml	W1001-4
	6 ml	W1001-6
	10 ml	W1001-10
<b>Zymo-Spin™ IIC Columns</b>	50	C1011-50
	250	C1011-250
<b>Zymo-Spin™ IICG Columns</b>	50	C1006-50-G
	250	C1006-250-G
<b>Collection Tubes</b>	50	C1001-50
	500	C1001-500
	1000	C1001-1000

**ZYMO RESEARCH CORP.**



**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple*

---

**ZYMO RESEARCH CORP.**

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ [info@zymoresearch.com](mailto:info@zymoresearch.com) ▪ [www.zymoresearch.com](http://www.zymoresearch.com)