

GenePurgeDirect® DNA/RNA Releasing Agent

FFPE and Tissue

Version: 1.0

Revision date: 31-07-2014

Description

GenePurgeDirect® is composed of proprietary polymeric materials that quickly facilitate the release of nucleic acids from cells in a form suitable for PCR. By segregating inhibitors that are released during lysis as well as any preservation agents that may interfere with amplification GenePurgeDirect® provides amplifiable nucleic acids from minute amounts of material. The protocols for Tissue/FFPE tissue are compiled from user-developed methods see the listed references for additional information.

Protocols

I.!!GenePurgeDirect®!Protocol!for!Fresh/Frozen!Tissue!

Reagents: Sterile DI water
1X TE
Proteinase K (optional)
GenePurgeDirect®

Part I: Tissue Homogenization

1. Cut a 1mm³ thick section from either fresh or frozen tissue that has been rinsed with sterile water to remove any surface contamination.
2. Place the section into the bottom of a 1.5ml tube.
3. Add 25µl of 1X TE to the tube containing the sectioned tissue.
4. Mince the section of tissue by pushing a pestle against the tissue and twisting the pestle to compress the tissue against the walls of the tube. Ten thrusts with the pestle should be sufficient.
* Optional: 10 units of Proteinase K can be added and the tissue digested for 113 hours at 55°C.
5. Place the tube containing the homogenate at 4°C if subsequent sections need to be homogenized.

Part II: GenePurgeDirect® Treatment

6. Flick the homogenate tube 3 times then transfer 1µl of the homogenate to a 0.5ml amplification tube that is PCR compatible with your thermal cycler.
7. Resuspend the GenePurgeDirect® mixture by vortexing 2-3 seconds or inverting 5-10 times.
8. Add 20µl of GenePurgeDirect® suspension to the 1µl of homogenate in the PCR tube and tightly close the tube lid. Pulse vortex briefly to mix.
9. Place samples onto thermal cycler with a heated lid with the following GenePurgeDirect® program:

<u>Step</u>	<u>Temperature</u>	<u>Time</u>
1.	65°C	30 sec.
2.	8°C	30 sec.
3.	65°C	90 sec.
4.	97°C	180 sec.
5.	8°C	60 sec.
6.	65°C	180 sec.
7.	97°C	60 sec.
8.	65°C	60 sec.
6.	80°C	hold

10. Once program is completed centrifuge sample tubes at 5 000xg for 1 minute.
11. Transfer the supernatant to a new tube for use as the template for PCR. Use 1110µl of supernatant per 201100µl amplification reactions. *PCR can be performed directly in the GenePurgeDirect® treatment tube; add amplification reagents for a final volume of 100µl.
12. Perform amplification reaction according to your optimized protocol.

II. GenePurgeDirect®!Protocol!for!Paraffin!Embedded!Tissue

Reagents: Sterile DI Water Xylene (Acetone if lipid rich tissue) Ethanol (100% 70% 30%)
1X TE
Proteinase K (optional)
GenePurgeDirect®

Part I: Deparaffinize tissue (either slide mounted or section form a block)

1. Perform 2 x 15 minute washed in xylene or equivalent by completely immersing the tissue section or slide in the wash agent. *Use 100% acetone wash 2 x 5 minutes for lipid rich tissue like brain tissue.
2. Perform the following sequential washes by completely immersing tissue section or slide in the wash agent:



Part II: GenePurgeDirect® Treatment

- Resuspend the GenePurgeDirect® mixture by vortexing 2-3 seconds or inverting 5-10 times.
- Combine tissue and 20ul of GenePurgeDirect® into a 1.5 ml tube. If tissue was slide mounted then use approximately 1mm² to 9mm² tissue scraping into the tube.
- Using a pestle grind the section of tissue by pushing a pestle against the tissue and twisting the pestle to compress the tissue against the walls of the tube. Ten thrusts with the pestle should be sufficient.
*Optional: 10 units of Proteinase K can be added and the tissue digested for 30 minutes at 55°C.
- Transfer as much of the tissue homogenate as possible to a new 0.5ml PCR tube.
- Proceed with either the thermal cycler lysis protocol (as described on page 1) or the microwave lysis protocol (below).
- After the lysis program is complete add 25µl of 1XTE to the treated tissue.
- Centrifuge sample tube at 5 000xg for 5 minutes.
- Transfer the supernatant to a fresh tube for use as a template for PCR. Use 1J10µl of supernatant per 20J100µl amplification reactions.
- Perform amplification reaction according to your optimized protocol.

Microwave!Lysis!Protocol:!

We have found that the microwave treatment of specimens affords a rapid sample preparation and facilitates the amplification of the more intractable types of specimens.

A.!!Evaluation!of!microwave!!

Perform the following experiment to determine the optimal conditions for your tubes and microwave.

- Place 40µl DI water in the same type of tube that you will be using for GenePurgeDirect® treatment.
- Overlay each tube with mineral oil to prevent evaporation.
- Close the tubes place in microwave safe rack (polyethylene or propylene) and heat on high for 5 minutes.
- If any caps pop or tubes distort in any manner then place a separate beaker in the microwave with 150ml of room temperature DI water and repeat the above 3 steps the beaker of water serves as a heat ballast.
- If tubes open or distort reduce the power by 10% increments and increase time by 1Jminute increments repeating step 4 until tubes no longer open or distort.

Note:!!Make!sure!the!racks!used!in!this!procedure!are!MICROWAVE!SAFE!!

B.!!Microwave!Protocol!

- Perform microwave procedure above for time and power conditions
- Place 1µl of specimen with 20µl of GenePurgeDirect® into either a 0.5ml PCR tube or 1.5ml tube.
- Vortex the tubes containing specimen and GenePurgeDirect® for ~10 seconds.
- Overlay with mineral oil to prevent samples from evaporating.
- Place the closed tubes in a microwave safe polyethylene or propylene rack. Make sure that the lids are loosely closed. If lids are closed too tightly tubes could rupture.
- Place the rack in a microwave oven and heat at maximum power setting(setting should be based on the microwave evaluation results) for 5J7 minutes. Typically 5 minutes if wattage is 900 or higher and 7 minutes if wattage is 500.
- Remove rack from microwave and centrifuge the tubes at 5000xg for 5 minutes. After centrifuging samples remove supernatant and use as DNA template.
- Perform the amplification reaction

References:!

Annis SL Dastoor FP Ziel H Daszak P Longcore JE. A DNAbased assay identifies Batrachochytrium dendrobatidis in amphibians. Journal of Wildlife Diseases. 2004 Jul;40(3):420J8.

Flore O Rafii S Ely S O'Leary JJ Hyjek EM Cesarman E. Transformation of primary human endothelial cells by Kaposi's sarcomaAssociated herpesvirus. Nature. 1998 Aug;394(6693):588J92.

Nakata B Wang YQ Yashiro M Nishioka N Tanaka H Ohira M Ishikawa T Nishino H Hirakawa K. Prognostic value of microsatellite instability in resectable pancreatic cancer. Clinical Cancer Research. 2002 Aug;8(8):2536J40.

Toledo LF Haddad CFB Carnaval ACOQ Britto FB. A Brazilian anuran (Hylodes magalhaesi : Leptodactylidae) infected by Batrachochytrium dendrobatidis : a conservation concern. Amphibian and Reptile Conservation 2006; 4(1):17J21.

India Contact

Biolinkk

Building No. 1333 Office No. 412 Durga Chamber D. B. Gupta Road, Karol Bagh, New Delhi-11005

Ph: : 011-45035753, 011-28753658 Email: customerservice@biolinkk.com