


NimaPOP™

High Quality Genetic Analyzer Polymers



Capillary Electrophoresis Polymers Designed for Applied Biosystems™ Genetic Analyzers


- > High Resolution & Performance
- > Long Shelf Life
- > Compatible with all Applied Biosystems™ Genetic Analyzers
- > No need for changes in any settings/calibrations

Introducing Nimapop™

NimaGen is proud to introduce two newly developed separation matrixes using state-of-the-art polymer technology: NimaPOP™4, NimaPOP™6 and NimaPOP™7.

These new polymers have an excellent dynamic coating and separating ability with all different applications. NimaPOP™ polymers can be used without any requirement for changes in run protocol, conditions or spectral calibrations.

Developed with the latest polymer chemistry technology, NimaPOP™ demonstrates an increased stability and resolution.



Compatibility

The polymers are available in different formats:

- Click-in bottles of 5, 10 and 28 ml for 3130 and 3730 series Genetic Analyzers
- Vials of 5 and 10 ml for 310 and 3100 series Genetic Analyzers
- Pouches for 384 and 960 samples for 3500 series Genetic Analyzers

NimaGen Polymers for Capillary Electrophoresis based Genetic Analyser

Analyser Applied BioSystem	NimaPOP™-4	NimaPOP™-6	NimaPOP™-7
NimaPOP™ for 3130 series, 5 ml	NIP4-005	NIP6-005	NIP7-005
NimaPOP™ for 3130 series, 10 ml	NIP4-010	NIP6-010	NIP7-010
NimaPOP™ for 310/3100 series, 5 ml	NIP4-105	NIP6-105	N/A
NimaPOP™ for 310/3100 series, 10 ml	NIP4-110	NIP6-110	N/A
NimaPOP™ for 3730 series, 28 ml	N/A	NIP6-028	NIP7-028
NimaPOP™ for 3730 series, 10 x 28 ml	N/A	NIP6-280	NIP7-280
NimaPOP™ for 3500 series, 384 samples	NIP4-384	NIP6-384	NIP7-384
NimaPOP™ for 3500 series, 960 samples	NIP4-960	NIP6-960	NIP7-960

Trace Score & QV20+ Reading Length Comparison NimaPOP™ with POP7

	POP7	NimaPOP™-7 Polymer
Score (ave.)	31.83	32.28
QV20+ Length (ave.)	590.0	601.4

Data is obtained by analyzing ~3,900 samples, including PCR products and plasmids

NimaPOP™ polymers are compatible for use with existing sequencing analysis software including KB™ basecaller, Phred and PeakTrace and deliver excellent performance with all BigDye®, BrighDye™ and BrilliantDye™ Cycle sequencing kits, as well as all Fragment Analysis applications.

Extended 3500 On-Instrument lifespan

The RFID labels supplied with the NimaPOP™ polymers and NimaPOP™ buffers have been formulated with an increased “on-instrument” lifespan, doubling its on instrument use when compared to the original product from 7 to 14 days.

NimaPOP™ for 3500 series Genetic Analyser

Polymers

Available in preformulated pouches with Radio Frequency Identification (RFID) labels, NimaPOP™ polymers for the 3500 series Genetic Analyser, that are compatible for direct connection to the instrument.

Running buffer for 4 complete refills

The NimaPOP™ 10x running buffer is available in a quantity of 60ml, enabling you to make up 600ml buffer.

The process is easy; just empty the used containers, rinse and refill with the new buffer.

Replace the used RFID label with the new, self-adhesive, label supplied with the NimaPOP™ buffer. With every bottle of 10x NimaPOP running buffer, you'll receive 4 new RFID labels for the ABC and 4 new labels for the CBC container.



NimaGen Capillary Electrophoresis Running Buffers

Product	Part-no.	Quantity
NimaPOP™ 10x Running Buffer (with EDTA) <i>for use with NimaPOP™ polymers</i>	NIB-025	25 ml. bottle
NimaPOP™ 10x Running Buffer (with EDTA) <i>for use with NimaPOP™ polymers</i>	NIB-100	100 ml. bottle
NimaPOP™ 10x Running Buffer (with EDTA) <i>for use with NimaPOP™ polymers</i>	NIB-500	500 ml. bottle
NimaPOP™ 10x Running Buffer for 3500 series <i>(Data Collection Software v1.0)</i> <i>incl. 4 x ABC & 4 x CBC self adhesive RFID chips</i>	NIB-3500	60 ml. bottle <i>10x buffer for 4 refills in total)</i>

GET
your **Quote**
NOW

Contact

Visiting address

NimaGen B.V.
Lagelandseweg 56
6545 CG Nijmegen
The Netherlands

India Contact

Biolinkk

**Building No 1333 Office No 412 ,Durga Chamber D.B.Gupta Road ,Karol Bagh,
New Delhi-110005**



Discover what we can do for you

www.nimagen.com