



Small, but convenient!

Mini Cell Strainers

For more information : http://www.funakoshi.co.jp/exports_contents/80727

Small and convenient cell strainers, which fit for 1.5 / 2.0 ml tubes, 5 ml sampling tubes and 15 ml conical tubes.

Mini Cell Strainer (for 5 ml sampling tubes and 15 ml conical tubes)



5 ml tube 15 ml tube 5ml tube



25 μ m 40 μ m 70 μ m



Features

- ◆ Facilitate speedy cell separation from primary tissues or cell clumps
- ◆ Slits on the rim avoid overflow and help speedy filtration
- ◆ Fit for 5 ml sampling and 15 ml conical tubes
- ◆ Color-coded according to mesh pore sizes (25 μ m : White, 40 μ m : Blue, 70 μ m : Clear)
- ◆ Materials : Frame ; PP, Mesh ; Nylon
- ◆ EOG sterilized and individually wrapped
- ◆ Made in Japan

Mini Cell Strainer II (for 1.5 / 2.0 ml or 5 ml ($\phi 12 \times 75$ mm) tubes)



1.5 ml tube



2 ml tube



5 ml tube
($\phi 12 \times 75$ mm)



Features

- ◆ For prefiltration of samples before flow cytometry or mass spectrometry
- ◆ Slits on the rim avoid overflow and help speedy filtration
- ◆ Fit for 1.5 / 2.0 ml tubes ($\phi 12$ mm) and 5 ml ($\phi 12 \times 75$ mm) sampling tubes
- ◆ EOG sterilized and individually wrapped
- ◆ Made in Japan

Product Information

■ Mini Cell Strainer

[Manufacturer : HIT]

Product Name	Size	Code
Mini Cell Strainer, 25 μ m, White	100 pieces	HT-AMS-12502
Mini Cell Strainer, 40 μ m, Blue	100 pieces	HT-AMS-14002
Mini Cell Strainer, 70 μ m, Clear	100 pieces	HT-AMS-17002

■ Mini Cell Strainer II

[Manufacturer : HIT]

Product Name	Size	Code
Mini Cell Strainer II, 40 μ m, Blue	100 pieces	HT-AMS-04002
Mini Cell Strainer II, 40 μ m, Blue with 5 ml ($\phi 12 \times 75$ mm) tube set (100 pieces each)	1 set	HT-AMS-04008

NOTE

* All products shown here are for research use only, not for diagnostic use.
* Products might be changed for improvement without notice.

* All brand and product names are registered as trademarks of their respective companies.
* Please contact your local distributors for orders, quote request and any inquiry.

Your Local Distributor

Biolinkk

Building NO: 1333, 4th Floor, Office No:412,
Durga Chamber, Desh Bandhu Gupta Road,
Karol Bagh, New Delhi- 110005, India
Phone: +91-11-45035753
email: customerservice@biolinkk.com



funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,
Tokyo 113-0033, JAPAN
Phone: +81-3-5684-6296
Fax: +81-3-5684-6297
Email: export@funakoshi.co.jp

HIT-1702-F02