


Convenient and Safe! Disposable Homogenizer BioMasher II

For more information : http://www.funakoshi.co.jp/exports_contents/46114

BioMasher II is a convenient disposable homogenizer. This is suitable for extraction of DNA, RNA and protein from animal, insect and plant tissues.




Features

- High homogenization efficiency
 - By texture on both pestle and inside of tube.
- High recovery
- Heat can be controlled
- No centrifuge
- Pestle has sample rebound structure.
- PP : Autoclavable / POM : non-Autoclavable
- EOG-Sterilized and non-sterilized product are available.
- Much higher homogenization efficiency by using PowerMasher II (Pestle Motor).

Handling procedure

1. Put samples and buffer into tube and insert pestle.
2. Grinding samples by pestle by hand or Power Masher II.



PowerMasher II

PowerMasher II is a pestle motor, optimized for BioMasher II, III, V and SP series.

Features

- Increase homogenization efficiency dramatically.
- Works on two AA batteries (Accessories).



Check a movie at
http://www.funakoshi.co.jp/exports_contents/80608

Product Information

■ BioMasher II

[Manufacturer : NIP]

| Product Name | Size | Code |
|-------------------------------|----------|--------|
| BioMasher II (non-sterilized) | 100 sets | 320102 |
| BioMasher II (EOG-sterilized) | 100 sets | 320103 |

■ PowerMasher II (Pestle motor)

[Manufacturer : NIP]

| Product Name | Size | Code |
|----------------|---------|--------|
| PowerMasher II | 1 piece | 891300 |

NOTE

- * All products on here are for research use only, not for diagnostic use.
- * Products might be changed for improvement without notice.

- * All brand names and product names are registered trademarks of their respective companies.

- * Please contact your local distributors for orders, quote request and any inquiry.

Your Local Distributor

Biolinkk

Building NO: 1333, 4th Floor, Office No:412,
 Durga Chamber, Desh Bandhu Gupta Road,
 Karol Bagh, New Delhi- 110005, India
 Phone: +91-11-45035753
 email: customerservice@biolinkk.com



funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,
 Tokyo 113-0033, JAPAN
 Phone: +81-3-5684-6296
 Fax: +81-3-5684-6297
 Email: export@funakoshi.co.jp