

04/09/2018

Kit Components

Product code	Description
--------------	-------------

D4080	Select-a-Size DNA Clean & Concentrator Kit (25 preps)
--------------	------------------------------------------------------------------

Components:

D4080-1-15	Select-a-Size DNA Binding Buffer
D4003-2-6	DNA Wash Buffer
D3004-4-1	DNA Elution Buffer



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

1 Identification

- **Product identifier**
- **Trade name:** *Select-a-Size DNA Binding Buffer*
- **Article number:** D4080-1-15, D4080-1-120
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:**
Integrated Sciences Pty.LTD.,
54 Gibbes Street, Chatswood NSW 2067,
Australia
Phone: 02 9417 7866
Fax: 02 9417 5066
tech@integratedsci.com.au
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS05 Corrosion

Skin Corr. 1C H314 Causes severe skin burns and eye damage.
 Eye Dam. 1 H318 Causes serious eye damage.



GHS07

Acute Tox. 4 H302 Harmful if swallowed.
 Acute Tox. 4 H312 Harmful in contact with skin.
 Acute Tox. 4 H332 Harmful if inhaled.

.....

Aquatic Chronic 3 H412 Harmful to aquatic life with long lasting effects.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS05, GHS07
- **Signal word** Danger

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 1)

· **Hazard-determining components of labeling:**

guanidinium thiocyanate

· **Hazard statements**

Harmful if swallowed, in contact with skin or if inhaled.

Causes severe skin burns and eye damage.

Harmful to aquatic life with long lasting effects.

· **Precautionary statements**

Do not breathe dusts or mists.

Wash thoroughly after handling.

Do not eat, drink or smoke when using this product.

Use only outdoors or in a well-ventilated area.

Avoid release to the environment.

Wear protective gloves/protective clothing/eye protection/face protection.

If swallowed: Call a poison center/doctor if you feel unwell.

If swallowed: Rinse mouth. Do NOT induce vomiting.

If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Immediately call a poison center/doctor.

Specific treatment (see on this label).

Take off contaminated clothing and wash it before reuse.

Wash contaminated clothing before reuse.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



Health = 3

Fire = 0

Reactivity = 0

· **HMIS-ratings (scale 0 - 4)**



Health = *3

Fire = 0

Reactivity = 0

· **Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Mixture of the substances listed below with nonhazardous additions.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 2)

Dangerous components:

Table with 2 columns: CAS: 593-84-0 guanidinium thiocyanate and ≤100%

4 First-aid measures

- Description of first aid measures
General information: Immediately remove any clothing soiled by the product.
After inhalation: Supply fresh air.
After skin contact: Immediately wash with water and soap.
After eye contact: Rinse opened eye for several minutes.
After swallowing: Rinse mouth.
Information for doctor: Most important symptoms and effects, both acute and delayed.
Indication of any immediate medical attention and special treatment needed.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: CO2, extinguishing powder or water spray.
Special hazards arising from the substance or mixture: Products of thermal decomposition.
Advice for firefighters
Protective equipment: Wear self-contained breathing apparatus.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures: Wear self-contained breathing apparatus.
Environmental precautions: Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 3)

- Methods and material for containment and cleaning up: Absorb with liquid-binding material... Use neutralizing agent... Dispose contaminated material... Ensure adequate ventilation.
Reference to other sections: See Section 7 for information on safe handling... See Section 8 for information on personal protection equipment... See Section 13 for disposal information.
Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number (593-84-0), chemical name (guanidinium thiocyanate), and concentration (0.98 mg/m3)

PAC-2:

Table with 3 columns: CAS number (593-84-0), chemical name (guanidinium thiocyanate), and concentration (11 mg/m3)

PAC-3:

Table with 3 columns: CAS number (593-84-0), chemical name (guanidinium thiocyanate), and concentration (65 mg/m3)

7 Handling and storage

- Handling:
Precautions for safe handling: Ensure good ventilation/exhaustion at the workplace. Prevent formation of aerosols.
Information about protection against explosions and fires: No special measures required.
Conditions for safe storage, including any incompatibilities: Store in cool, dry place. Store in well-ventilated location.
Storage:
Requirements to be met by storerooms and receptacles: No special requirements.
Information about storage in one common storage facility: Do not store together with acids or strong oxidizers.
Further information about storage conditions: Keep receptacle tightly sealed.
Specific end use(s): Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters: Work under a chemical fume hood when using this product. Ensure eyewash station and safety showers are readily accessible.
Components with limit values that require monitoring at the workplace: The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
Additional information: The lists that were valid during the creation were used as basis.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 4)

· **Exposure controls**

The appropriate protective equipment under anticipated circumstances of use include lab-coat, safety glasses with side-shields and gloves.

· **Personal protective equipment:**

· **General protective and hygienic measures:**

- Keep away from foodstuffs, beverages and feed.
- Immediately remove all soiled and contaminated clothing.
- Wash hands before breaks and at the end of work.
- Avoid contact with the eyes and skin.

· **Breathing equipment:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:**



Tightly sealed goggles

9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· Form:	Liquid
· Color:	Clear
· Odor:	No information available
· Odor threshold:	Not determined.

· **pH-value:** Not determined.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 5)

Table with 2 columns: Property and Value. Rows include: Change in condition (Melting point, Boiling point), Flash point, Flammability, Ignition temperature, Decomposition temperature, Auto igniting, Danger of explosion, Explosion limits, Vapor pressure, Density, Relative density, Vapor density, Evaporation rate, Solubility in / Miscibility with Water, Partition coefficient, Viscosity, Solvent content, VOC content, Solids content, and Other information.

10 Stability and reactivity

- Reactivity: No further relevant information available.
Chemical stability: This product is normally stable under anticipated circumstances of use and storage.
Thermal decomposition / conditions to be avoided: Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon nitrogen and sulfur.
Possibility of hazardous reactions: No dangerous reactions known.
Conditions to avoid: Avoid exposing product to extreme temperatures or incompatible chemicals.
Incompatible materials: Acids and strong oxidizers.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 6)

- Hazardous decomposition products:
Product will not undergo self-decomposition, so no such products will be generated.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:

LD/LC50 values that are relevant for classification:

CAS: 593-84-0 guanidinium thiocyanate

Oral LD50 593 mg/kg (rat)

Primary irritant effect:

on the skin: Caustic effect on skin and mucous membranes.

on the eye: Strong caustic effect.

Sensitization: No sensitizing effects known.

Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations:

Harmful

Corrosive

Swallowing will lead to a strong caustic effect on mouth and throat and to the danger of perforation of esophagus and stomach.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

Toxicity

Aquatic toxicity:

CAS: 593-84-0 guanidinium thiocyanate

EC50 42.4 mg/kg (daphnia)

Persistence and degradability No further relevant information available.

Behavior in environmental systems:

Bioaccumulative potential No further relevant information available.

Mobility in soil No further relevant information available.

Additional ecological information:

General notes:

Water hazard class 2 (Self-assessment): hazardous for water

Do not allow product to reach ground water, water course or sewage system.

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 7)

Must not reach bodies of water or drainage ditch undiluted or unneutralized. Danger to drinking water if even small quantities leak into the ground.

Results of PBT and vPvB assessment

- PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects: No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Dispose of contents in accordance with local/regional/national, and international recommendations.
Uncleaned packagings:
Recommendation: Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Description and Value. Rows include UN-Number (UN1760), UN proper shipping name (Corrosive liquids, n.o.s. (guanidinium thiocyanate)), Transport hazard class(es) (8 Corrosive substances), and Packing group (III).

(Contd. on page 9)



Safety Data Sheet acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 8)

Table with hazard information: Environmental hazards, Special precautions for user, Danger code, EMS Number, Stowage Category, Stowage Code, Transport in bulk, Transport/Additional information, DOT, Quantity limitations, IMDG, Limited quantities, Excepted quantities, UN Model Regulation.

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture Sara

Table with regulatory information: Section 355, Section 313, TSCA, Proposition 65, Chemicals known to cause cancer, Chemicals known to cause reproductive toxicity for females, Chemicals known to cause reproductive toxicity for males.

(Contd. on page 10)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 9)

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

· Hazard pictograms GHS05, GHS07

· Signal word Danger

· Hazard-determining components of labeling:

guanidinium thiocyanate

· Hazard statements

Harmful if swallowed, in contact with skin or if inhaled.

Causes severe skin burns and eye damage.

Harmful to aquatic life with long lasting effects.

· Precautionary statements

Do not breathe dusts or mists.

Wash thoroughly after handling.

Do not eat, drink or smoke when using this product.

Use only outdoors or in a well-ventilated area.

Avoid release to the environment.

Wear protective gloves/protective clothing/eye protection/face protection.

If swallowed: Call a poison center/doctor if you feel unwell.

If swallowed: Rinse mouth. Do NOT induce vomiting.

If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Immediately call a poison center/doctor.

Specific treatment (see on this label).

Take off contaminated clothing and wash it before reuse.

Wash contaminated clothing before reuse.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

(Contd. on page 11)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 11/11

Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 03/20/2017

Trade name: Select-a-Size DNA Binding Buffer

(Contd. of page 10)

· **Department issuing SDS:**

Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614
USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 04/09/2018 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Acute Tox. 4: Acute toxicity – Category 4

Skin Corr. 1C: Skin corrosion/irritation – Category 1C

Eye Dam. 1: Serious eye damage/eye irritation – Category 1

Aquatic Chronic 3: Hazardous to the aquatic environment - long-term aquatic hazard – Category 3

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/16/2016

1 Identification

- **Product identifier**
- **Trade name:** DNA Wash Buffer
- **Article number:** D4003-2-6, D4003-2-24, D4003-2-48
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:**
Integrated Sciences Pty.LTD.,
54 Gibbes Street, Chatswood NSW 2067,
Australia
Phone: 02 9417 7866
Fax: 02 9417 5066
tech@integratedsci.com.au
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer

vPvB: Not applicable.

(Contd. of page 1)

3 Composition/information on ingredients

- Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.
Dangerous components: Void

4 First-aid measures

- Description of first aid measures
General information: No special measures required.
After inhalation: Supply fresh air; consult doctor in case of complaints.
After skin contact: Generally the product does not irritate the skin.
After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.
Indication of any immediate medical attention and special treatment needed
No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam. Use fire fighting measures that suit the environment.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures Wear protective clothing.
Environmental precautions: Dilute with plenty of water.
Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Reference to other sections
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
Protective Action Criteria for Chemicals

PAC-1:

CAS: 1185-53-1 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride

12 mg/m³

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer

(Contd. of page 2)

CAS: 6381-92-6	Edetate Disodium, Dihydrate	30 mg/m ³
· PAC-2:		
CAS: 1185-53-1	2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride	130 mg/m ³
CAS: 6381-92-6	Edetate Disodium, Dihydrate	330 mg/m ³
· PAC-3:		
CAS: 1185-53-1	2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride	790 mg/m ³
CAS: 6381-92-6	Edetate Disodium, Dihydrate	2,000 mg/m ³

7 Handling and storage

- **Handling:**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Not required.
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer

(Contd. of page 3)

- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Goggles recommended during refilling.

9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

- **Form:** Liquid
- **Color:** Clear
- **Odor:** Odorless
- **Odor threshold:** Not determined.

· **pH-value:** Not determined.

· **Change in condition**

- **Melting point/Melting range:** Undetermined.
- **Boiling point/Boiling range:** Undetermined.

· **Flash point:** Not applicable.

· **Flammability (solid, gaseous):** Not applicable.

· **Ignition temperature:**

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Product is not selfigniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

- **Lower:** Not Applicable
- **Upper:** Not Applicable

· **Vapor pressure:** Not determined.

- **Density:** Not determined.
- **Relative density** Not determined.
- **Vapor density** Not determined.
- **Evaporation rate** Not determined.

· **Solubility in / Miscibility with**

· **Water:** Fully miscible.

· **Partition coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

- **Dynamic:** Not determined.
- **Kinematic:** Not determined.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer

(Contd. of page 4)

Table with 2 columns: Property (Solvent content, VOC content, Solids content, Other information) and Value (0.00 %, 0.0 g/l / 0.00 lb/gl, 1.0 %, No further relevant information available).

10 Stability and reactivity

- Reactivity, Chemical stability, Thermal decomposition / conditions to be avoided, Possibility of hazardous reactions, Conditions to avoid, Incompatible materials, Hazardous decomposition products.

11 Toxicological information

- Information on toxicological effects, Acute toxicity, Primary irritant effect (on the skin, on the eye), Sensitization, Additional toxicological information.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity, Aquatic toxicity, Persistence and degradability.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer

(Contd. of page 5)

- Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes: Generally not hazardous for water
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Smaller quantities can be disposed of with household waste.
Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Description and Regulation status. Rows include UN-Number, UN proper shipping name, Transport hazard class(es), Packing group, Environmental hazards, Special precautions for user, Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code, and UN Model Regulation.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer

(Contd. of page 6)

15 Regulatory information

- *Safety, health and environmental regulations/legislation specific for the substance or mixture*
- *Sara*

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act):**

CAS: 1185-53-1 | 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride

· **TSCA new (21st Century Act) (Substances not listed)**

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** Void

· **Hazard pictograms** Void

· **Signal word** Void

· **Hazard statements** Void

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer

(Contd. of page 7)

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 04/09/2018 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/03/2015

1 Identification

- **Product identifier**
- **Trade name:** DNA Elution Buffer
- **Article number:** D3004-4-1, D3004-4-4, D3004-4-10, D3004-4-16, D3004-4-50
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 1)

3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of the substances listed below with nonhazardous additions.
- **Dangerous components:** Void

4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire fighting measures that suit the environment.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:**
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· **PAC-1:**

CAS: 1185-53-1 | 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride

12 mg/m³

(Contd. on page 3)



Safety Data Sheet

acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 2)

CAS: 6381-92-6	Edetate Disodium, Dihydrate	30 mg/m ³
· PAC-2:		
CAS: 1185-53-1	2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride	130 mg/m ³
CAS: 6381-92-6	Edetate Disodium, Dihydrate	330 mg/m ³
· PAC-3:		
CAS: 1185-53-1	2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride	790 mg/m ³
CAS: 6381-92-6	Edetate Disodium, Dihydrate	2,000 mg/m ³

7 Handling and storage

- **Handling:**
- **Precautions for safe handling:** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s):** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Not required.
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 3)

- Penetration time of glove material
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

- Form: Liquid
Color: Clear
Odor: Odorless
Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

- Melting point/Melting range: Undetermined.
Boiling point/Boiling range: Undetermined.

Flash point: Not applicable.

Flammability (solid, gaseous): Not applicable.

Ignition temperature:

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

- Lower: Not Applicable
Upper: Not Applicable

Vapor pressure: Not determined.

- Density: Not determined.
Relative density: Not determined.
Vapor density: Not determined.
Evaporation rate: Not determined.

Solubility in / Miscibility with

Water: Fully miscible.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

- Dynamic: Not determined.
Kinematic: Not determined.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 4)

· Solvent content:	
· VOC content:	0.00 % 0.0 g/l / 0.00 lb/gl
· Solids content:	2.0 %
· Other information	No further relevant information available.

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**
The product is not subject to classification according to internally approved calculation methods for preparations: When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

· **Carcinogenic categories**

· **IARC (International Agency for Research on Cancer)**

None of the ingredients is listed.

· **NTP (National Toxicology Program)**

None of the ingredients is listed.

· **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.

(Contd. on page 6)



Safety Data Sheet acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 5)

- Behavior in environmental systems: Bioaccumulative potential No further relevant information available. Mobility in soil No further relevant information available. Additional ecological information: General notes: Water hazard class 1 (Self-assessment): slightly hazardous for water Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system. Results of PBT and vPvB assessment: PBT: Not applicable. vPvB: Not applicable. Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods Recommendation: Smaller quantities can be disposed of with household waste. Uncleaned packagings: Recommendation: Dispose of container in accordance with local/regional/national and international recommendations. Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Hazard/Requirement and Status. Rows include UN-Number, UN proper shipping name, Transport hazard class(es), Packing group, Environmental hazards, Special precautions for user, Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code, and UN Model Regulation.



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 6)

15 Regulatory information

- Safety, health and environmental regulations/legislation specific for the substance or mixture
- Sara

· **Section 355 (extremely hazardous substances):**
None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**
None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act):**
CAS: 1185-53-1 | 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride

- TSCA new (21st Century Act) (Substances not listed)
- Proposition 65

· **Chemicals known to cause cancer:**
None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**
None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**
None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**
None of the ingredients is listed.

- Carcinogenic categories

· **EPA (Environmental Protection Agency)**
None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**
None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**
None of the ingredients is listed.

- GHS label elements Void
- Hazard pictograms Void
- Signal word Void
- Hazard statements Void
- Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Department issuing SDS:**
Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet
acc. to OSHA HCS

Printing date 04/09/2018

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 7)

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 04/09/2018 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

US