


# DNA Clean & Concentrator® Kits


- **Fast & Simple:** Clean and concentrate DNA in 2 minutes.
- **Highly Concentrated:** Elute in as little as 6 µl.
- **Ultra-Pure:** Ready for Next-Gen Sequencing, PCR, ligations, etc.

## Simple, Innovative Workflow




Product	Cat. No.	Size
DNA Clean & Concentrator® -5 (uncapped columns)	D4003 D4003T D4004	50 preps 10 preps 200 preps
DNA Clean & Concentrator® -25 (uncapped columns)	D4005 D4006	50 preps 200 preps
ZR-96 DNA Clean & Concentrator®-5	D4023 D4024	2 x 96 preps 4 x 96 preps
Genomic DNA Clean & Concentrator®-10	D4010 D4011	25 preps 100 preps
Genomic DNA Clean & Concentrator®-25	D4064 D4065	25 preps 100 preps
ZR-96 Genomic DNA Clean & Concentrator®-5	D4066 D4067	2 x 96 preps 4 x 96 preps

## Recover Ultra-Pure, High-Quality DNA



High efficiency DNA recovery using the DNA Clean & Concentrator®-5 (DCC®-5). Equivalent amounts of DNA were purified using the DCC®-5 or MinElute® Kit (Supplier Q) and analyzed by electrophoresis in a 0.8% (w/v) agarose/TAE/EtBr gel.

## Improved Recovery of Large DNA



High molecular weight DNA is efficiently purified using the Genomic DNA Clean & Concentrator®-10. Porcine gDNA (~35-50 kb), T4 phage DNA (170 kb), and lambda (λ) phage DNA (48.5 kb) were purified (in duplicate) from input material using the Genomic DCC®. Eluted DNAs were analyzed in a 0.8% (w/v) TAE/agarose/EtBr gel. The size marker "M" is a 1 kb ladder (Zymo Research).

