

# DNA/RNA Clean-Up and Size Selection


## DNA/RNA Clean & Concentrator™

(library preparations, PCR, restriction digests and other enzymatic reactions, etc.)


- **Tunable:** Size selection can be tuned from 100 bp to 1000 bp with left, right, or double size selection.
- **Ultra-Pure:** 10 µl elutions are ready for Next-Gen Sequencing, etc.
- **Automation Ready:** Scripts and automation support readily available.

### Easy Size Selection

Select-a-Size DNA Clean & Concentrator® MagBead Kit allows for adjustable size selection. Exemplary size selections (using 2 µg of sonicated DNA) were analyzed using the Agilent 2200 TapeStation® system.




### High DNA Yield



Select-a-Size DNA Clean & Concentrator® MagBead shows efficient recovery at different concentrations. DNA recovery ( $\geq 80\%$ ) is consistent from low to high amounts of genomic DNA input (n=3).

### Small RNA Separation



Small RNA (<200 nt) is separated from large RNA (>200 nt) using the RNA Clean & Concentrator™ MagBead. Size was analyzed using the Agilent 2200 TapeStation® system.

Product	Cat. No.	Size
Select-a-Size DNA Clean & Concentrator® MagBead Kit	D4084 D4085	10 ml 50 ml
RNA Clean & Concentrator™ MagBead Kit	R1081	10 ml

[info@zymoresearch.com](mailto:info@zymoresearch.com)

[www.zymoresearch.com](http://www.zymoresearch.com)

Toll Free: (888) 882-9682