

08/22/2017

Kit Components

Product code	Description
--------------	-------------

D4306	ZymoBIOMICS 96 MagBead DNA Kit (2 x 96 preps)
--------------	--

Components:

D4302-2-250	ZymoBIOMICS MagBinding Buffer
D4302-3-100	ZymoBIOMICS MagWash1
D4302-4-200	ZymoBIOMICS MagWash 2
D4302-5-1, D4302-5-10	ZymoBIOMICS DNase/RNase Free Water
D4302-6-12	ZymoBIOMICS™ MagBinding Beads



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS MagBinding Buffer
- **Article number:** D4302-2-250
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS05 Corrosion

Skin Corr. 1C H314 Causes severe skin burns and eye damage.

Eye Dam. 1 H318 Causes serious eye damage.



GHS07

Acute Tox. 4 H302 Harmful if swallowed.

Aquatic Chronic 3 H412 Harmful to aquatic life with long lasting effects.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS05, GHS07
- **Signal word** Danger
- **Hazard-determining components of labeling:**
guanidinium thiocyanate
- **Hazard statements**
Harmful if swallowed.
Causes severe skin burns and eye damage.
Harmful to aquatic life with long lasting effects.
- **Precautionary statements**
Do not breathe mist/vapours/spray.
Wear protective gloves/protective clothing/eye protection/face protection.
Avoid release to the environment.

(Contd. on page 2)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 1)

Wash thoroughly after handling. Do not eat, drink or smoke when using this product. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Specific treatment (see on this label). IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Wash contaminated clothing before reuse. If swallowed: Rinse mouth. Do NOT induce vomiting. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

Classification system:

NFPA ratings (scale 0 - 4)



HMIS-ratings (scale 0 - 4)



Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

3 Composition/information on ingredients

Chemical characterization: Mixtures

Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: CAS number and component name, and a percentage column. Rows include guanidinium thiocyanate (≤50%) and glycerol (≤50%).

4 First-aid measures

Description of first aid measures

General information:

Immediately remove any clothing soiled by the product. Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 2)

- **After inhalation:**
Supply fresh air. If required, provide artificial respiration if trained to do so. Keep patient warm. Consult doctor if symptoms persist.
In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:**
Rinse mouth
DO NOT induce vomiting.
- **Information for doctor:**
· **Most important symptoms and effects, both acute and delayed** No further relevant information available.
· **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture**
Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon, nitrogen and sulfur.
- **Advice for firefighters**
- **Protective equipment:** Wear self-contained breathing apparatus for fighting fires involving this material

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**
Wear self-contained breathing apparatus for responding to non-incident release of this material in which there is the potential for inhalation of vapors, mists or sprays
Wear protective equipment. Keep unprotected persons away.
- **Environmental precautions:**
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Use neutralizing agent.
Dispose contaminated material as waste according to item 13.
Ensure adequate ventilation.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

(Contd. on page 4)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 3)

Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium thiocyanate (0.98 mg/m3) and glycerol (45 mg/m3).

PAC-2:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium thiocyanate (11 mg/m3) and glycerol (180 mg/m3).

PAC-3:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include guanidinium thiocyanate (65 mg/m3) and glycerol (1,100 mg/m3).

7 Handling and storage

- Handling:
Precautions for safe handling: Ensure good ventilation/exhaustion at the workplace.
Information about protection against explosions and fires: No special measures required.
Conditions for safe storage, including any incompatibilities: Store in cool, dry place.
Storage:
Requirements to be met by storerooms and receptacles: No special requirements.
Information about storage in one common storage facility: Do not store together with acids or strong oxidizers.
Further information about storage conditions: Keep receptacle tightly sealed.
Specific end use(s): Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters: Work under a chemical fume hood when using this product.
Components with limit values that require monitoring at the workplace: The following constituent is the only constituent of the product which has a PEL, TLV or other recommended exposure limit.

Table with 2 columns: Exposure limit type and value. Row 1: CAS: 56-81-5 glycerol. Row 2: PEL Long-term value: 15* 5** mg/m3 mist; *total dust **respirable fraction. Row 3: TLV TLV withdrawn-insufficient data human occup. exp.

Additional information: The lists that were valid during the creation were used as basis.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 4)

· Exposure controls

The appropriate protective equipment under anticipated circumstances of use include lab-coat, safety glasses with side-shields and gloves.

· Personal protective equipment:

· General protective and hygienic measures:

- Keep away from foodstuffs, beverages and feed.
- Immediately remove all soiled and contaminated clothing.
- Wash hands before breaks and at the end of work.
- Avoid contact with the eyes and skin.

· Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection:



Tightly sealed goggles

9 Physical and chemical properties

· Information on basic physical and chemical properties

· General Information

· Appearance:

- Form:** Liquid
- Color:** Clear
- Odor:** Mild
- Odor threshold:** Not determined.

· pH-value: Not determined.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 5)

Table with 2 columns: Property and Value. Rows include: Change in condition (Melting point, Boiling point), Flash point, Flammability, Ignition temperature, Decomposition temperature, Auto igniting, Danger of explosion, Explosion limits (Lower, Upper), Vapor pressure at 20 °C, Density, Relative density, Vapor density, Evaporation rate, Solubility in / Miscibility with Water, Partition coefficient (n-octanol/water), Viscosity (Dynamic, Kinematic), Solvent content (Organic solvents, VOC content), Solids content, and Other information.

10 Stability and reactivity

- Reactivity: No further relevant information available.
Chemical stability: This product is normally stable under anticipated circumstances of use and storage.
Thermal decomposition / conditions to be avoided: Products of thermal decomposition of this material would include hydrogen cyanide, ammonia, and oxides of carbon nitrogen and sulfur.
Possibility of hazardous reactions: No dangerous reactions known.
Conditions to avoid: Avoid exposing product to extreme temperatures or incompatible chemicals.
Incompatible materials: Acids and strong oxidizers.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 6)

- Hazardous decomposition products:
Product will not undergo self-decomposition, so no such products will be generated.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:

LD/LC50 values that are relevant for classification:

CAS: 593-84-0 guanidinium thiocyanate

Oral LD50 593 mg/kg (rat)

Primary irritant effect:

on the skin: Caustic effect on skin and mucous membranes.

on the eye: Strong caustic effect.

Sensitization: No sensitizing effects known.

Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations:

Harmful

Corrosive

Swallowing will lead to a strong caustic effect on mouth and throat and to the danger of perforation of esophagus and stomach.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

Toxicity

Aquatic toxicity:

CAS: 593-84-0 guanidinium thiocyanate

EC50 42.4 mg/kg (daphnia)

Persistence and degradability No further relevant information available.

Behavior in environmental systems:

Bioaccumulative potential No further relevant information available.

Mobility in soil No further relevant information available.

Additional ecological information:

General notes:

Water hazard class 2 (Self-assessment): hazardous for water

Do not allow product to reach ground water, water course or sewage system.

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 7)

Must not reach bodies of water or drainage ditch undiluted or unneutralized. Danger to drinking water if even small quantities leak into the ground.

Results of PBT and vPvB assessment

- PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects: No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Dispose of contents in accordance with local/regional/national, and international recommendations.
Uncleaned packagings:
Recommendation: Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Description and Value. Rows include UN-Number (UN1760), UN proper shipping name (CORROSIVE LIQUID, N.O.S. (guanidinium thiocyanate)), Transport hazard class(es) (8 Corrosive substances), and Packing group (III).

(Contd. on page 9)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 8)

Table with 2 columns: Hazard/Precaution and Description. Rows include Environmental hazards, Special precautions for user, Transport in bulk, and UN Model Regulation.

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Table with 2 columns: Regulation and Status. Rows include Section 355, Section 313, TSCA, Proposition 65, and various chemical listings.

(Contd. on page 10)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 9)

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms GHS05, GHS07

Signal word Danger

Hazard-determining components of labeling:

guanidinium thiocyanate

Hazard statements

Harmful if swallowed.

Causes severe skin burns and eye damage.

Harmful to aquatic life with long lasting effects.

Precautionary statements

Do not breathe mist/vapours/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

Avoid release to the environment.

Wash thoroughly after handling.

Do not eat, drink or smoke when using this product.

If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Immediately call a POISON CENTER/doctor.

Specific treatment (see on this label).

IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash contaminated clothing before reuse.

If swallowed: Rinse mouth. Do NOT induce vomiting.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Department issuing SDS:

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

(Contd. on page 11)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 11/11

Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagBinding Buffer

(Contd. of page 10)

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 08/22/2017 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Acute Tox. 4: Acute toxicity – Category 4

Skin Corr. 1C: Skin corrosion/irritation – Category 1C

Eye Dam. 1: Serious eye damage/eye irritation – Category 1

Aquatic Chronic 3: Hazardous to the aquatic environment - long-term aquatic hazard – Category 3

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS MagWash1
- **Article number:** D4302-3-100
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS02 Flame

Flam. Liq. 2 H225 Highly flammable liquid and vapor.



GHS07

Acute Tox. 4 H302 Harmful if swallowed.
Skin Irrit. 2 H315 Causes skin irritation.
Eye Irrit. 2A H319 Causes serious eye irritation.
STOT SE 3 H336 May cause drowsiness or dizziness.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS02, GHS07
- **Signal word** Danger
- **Hazard-determining components of labeling:**
guanidinium chloride
propan-2-ol
- **Hazard statements**
Highly flammable liquid and vapor.
Harmful if swallowed.
Causes skin irritation.
Causes serious eye irritation.
May cause drowsiness or dizziness.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 1)

· Precautionary statements

- Keep away from heat/sparks/open flames/hot surfaces. No smoking.
- Use explosion-proof electrical/ventilating/lighting/equipment.
- Avoid breathing dust/fume/gas/mist/vapors/spray
- Wear protective gloves / eye protection / face protection.
- Wear protective gloves / eye protection / face protection.
- Ground/bond container and receiving equipment.
- Use only non-sparking tools.
- Take precautionary measures against static discharge.
- Wash thoroughly after handling.
- Do not eat, drink or smoke when using this product.
- Use only outdoors or in a well-ventilated area.
- If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
- If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- Specific treatment (see on this label).
- IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
- IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- If skin irritation occurs: Get medical advice/attention.
- If eye irritation persists: Get medical advice/attention.
- Rinse mouth.
- In case of fire: Use for extinction: CO2, powder or water spray.
- Take off contaminated clothing and wash it before reuse.
- Store in a well-ventilated place. Keep container tightly closed.
- Store in a well-ventilated place. Keep cool.
- Store locked up.
- Dispose of contents/container in accordance with local/regional/national/international regulations.

· Classification system:

· NFPA ratings (scale 0 - 4)



· HMIS-ratings (scale 0 - 4)



· Other hazards

· Results of PBT and vPvB assessment

- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 2)

3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of the substances listed below with nonhazardous additions.

· **Dangerous components:**

CAS: 67-63-0	propan-2-ol	≤50%
CAS: 50-01-1	guanidinium chloride	≤50%

4 First-aid measures

- **Description of first aid measures**
- **General information:**
Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.
- **After inhalation:** In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** CO2, sand, extinguishing powder. Do not use water.
- **For safety reasons unsuitable extinguishing agents:** Water with full jet
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**
Wear protective equipment. Keep unprotected persons away.
- **Environmental precautions:**
Do not allow product to reach sewage system or any water course.
Prevent seepage into sewage system, workpits and cellars.
Inform respective authorities in case of seepage into water course or sewage system.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Dispose contaminated material as waste according to item 13.
Ensure adequate ventilation.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 3)

- Do not flush with water or aqueous cleansing agents
Reference to other sections
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include propan-2-ol (400 ppm) and guanidinium chloride (1.4 mg/m3).

PAC-2:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include propan-2-ol (2000* ppm) and guanidinium chloride (16 mg/m3).

PAC-3:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include propan-2-ol (12000** ppm) and guanidinium chloride (94 mg/m3).

7 Handling and storage

- Handling:
Precautions for safe handling
Information about protection against explosions and fires:
Conditions for safe storage, including any incompatibilities
Storage:
Requirements to be met by storerooms and receptacles:
Information about storage in one common storage facility:
Further information about storage conditions:
Specific end use(s) Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems:
Control parameters
Components with limit values that require monitoring at the workplace:
The following constituent is the only constituent of the product which has a PEL, TLV or other recommended exposure limit.
At this time, the remaining constituent has no known exposure limits.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 4)

CAS: 67-63-0 propan-2-ol

Table with 2 columns: PEL, REL, TLV and their corresponding values for propan-2-ol.

Ingredients with biological limit values:

CAS: 67-63-0 propan-2-ol

Table with 2 columns: BEI and its corresponding values and parameters for propan-2-ol.

Additional information: The lists that were valid during the creation were used as basis.

Exposure controls

Personal protective equipment:

General protective and hygienic measures:

- Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.

Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 5)

Eye protection:



Tightly sealed goggles

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

- Form: Liquid
Color: Clear
Odor: Alcohol-like
Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

- Melting point/Melting range: Undetermined.
Boiling point/Boiling range: Undetermined.

Flash point: 13 °C (55 °F)

Flammability (solid, gaseous): Not applicable.

Ignition temperature: 425 °C (797 °F)

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product is not explosive. However, formation of explosive air/vapor mixtures are possible.

Explosion limits:

- Lower: 2.0 Vol %
Upper: 12.0 Vol %

Vapor pressure at 20 °C (68 °F): 43 hPa (32 mm Hg)

- Density: Not determined.
Relative density: Not determined.
Vapor density: Not determined.
Evaporation rate: Not determined.

Solubility in / Miscibility with

Water: Not miscible or difficult to mix.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

Dynamic: Not determined.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 6)

Table with 2 columns: Property (Kinematic, Solvent content, Organic solvents, VOC content, Other information) and Value (Not determined, 49.0%, 49.0%, 490.0 g/l / 4.09 lb/gl, No further relevant information available).

10 Stability and reactivity

- Reactivity: No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions: No dangerous reactions known.
Conditions to avoid: No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:

LD/LC50 values that are relevant for classification:

CAS: 50-01-1 guanidinium chloride

Oral LD50 475 mg/kg (rat)

- Primary irritant effect:
on the skin: Irritant to skin and mucous membranes.
on the eye: Irritating effect.
Sensitization: No sensitizing effects known.

Additional toxicological information:

The product shows the following dangers according to internally approved calculation methods for preparations: Harmful, Irritant

Carcinogenic categories

IARC (International Agency for Research on Cancer)

CAS: 67-63-0 propan-2-ol

3

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

(Contd. on page 8)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 7)

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes: Not known to be hazardous to water.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation:
Dispose of contents in accordance with local/regional/national, and international recommendations.
Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.

14 Transport information

- UN-Number
DOT, IMDG, IATA UN1219
UN proper shipping name
DOT Isopropanol mixture
IMDG, IATA ISOPROPANOL (ISOPROPYL ALCOHOL) mixture
Transport hazard class(es)
DOT
Class 3 Flammable liquids



(Contd. on page 9)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 8)

Table with hazard information including Label (3), IMDG/IATA (Flammable liquids 3), Packing group (II), Environmental hazards (Not applicable), Special precautions (Warning: Flammable liquids), Transport in bulk (Not applicable), and UN Model Regulation (UN 1219 ISOPROPANOL MIXTURE, 3, II).

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture Sara

Section 355 (extremely hazardous substances):

None of the ingredients is listed.

Section 313 (Specific toxic chemical listings):

CAS: 67-63-0 | propan-2-ol

TSCA (Toxic Substances Control Act):

All ingredients are listed.

(Contd. on page 10)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 9)

Proposition 65

Chemicals known to cause cancer:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

CAS: 67-63-0 propan-2-ol

A4

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms GHS02, GHS07

Signal word Danger

Hazard-determining components of labeling:

guanidinium chloride

propan-2-ol

Hazard statements

Highly flammable liquid and vapor.

Harmful if swallowed.

Causes skin irritation.

Causes serious eye irritation.

May cause drowsiness or dizziness.

Precautionary statements

Keep away from heat/sparks/open flames/hot surfaces. No smoking.

Use explosion-proof electrical/ventilating/lighting/equipment.

Avoid breathing dust/fume/gas/mist/vapors/spray

Wear protective gloves / eye protection / face protection.

Wear protective gloves / eye protection / face protection.

Ground/bond container and receiving equipment.

Use only non-sparking tools.

Take precautionary measures against static discharge.

Wash thoroughly after handling.

Do not eat, drink or smoke when using this product.

Use only outdoors or in a well-ventilated area.

If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

(Contd. on page 11)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 11/12

Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 10)

Specific treatment (see on this label).

IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

If skin irritation occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

Rinse mouth.

In case of fire: Use for extinction: CO₂, powder or water spray.

Take off contaminated clothing and wash it before reuse.

Store in a well-ventilated place. Keep container tightly closed.

Store in a well-ventilated place. Keep cool.

Store locked up.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 08/22/2017 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

BEI: Biological Exposure Limit

Flam. Liq. 2: Flammable liquids – Category 2

Acute Tox. 4: Acute toxicity – Category 4

(Contd. on page 12)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 12/12

Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash1

(Contd. of page 11)

Skin Irrit. 2: Skin corrosion/irritation – Category 2
Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A
STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS MagWash 2
- **Article number:** D4302-4-200
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS02 Flame

Flam. Liq. 2 H225 Highly flammable liquid and vapor.



GHS07

Eye Irrit. 2A H319 Causes serious eye irritation.

STOT SE 3 H336 May cause drowsiness or dizziness.

- **Label elements**

- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

- **Hazard pictograms** GHS02, GHS07

- **Signal word** Danger

- **Hazard-determining components of labeling:**

propan-2-ol

ethanol

- **Hazard statements**

Highly flammable liquid and vapor.

Causes serious eye irritation.

May cause drowsiness or dizziness.

- **Precautionary statements**

Keep away from heat/sparks/open flames/hot surfaces. No smoking.

Use explosion-proof electrical/ventilating/lighting/equipment.

Avoid breathing dust/fume/gas/mist/vapors/spray

(Contd. on page 2)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 1)

Wear protective gloves / eye protection / face protection. Ground/bond container and receiving equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Wash thoroughly after handling. Use only outdoors or in a well-ventilated area. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/doctor if you feel unwell. If eye irritation persists: Get medical advice/attention. In case of fire: Use for extinction: CO2, powder or water spray. Store in a well-ventilated place. Keep container tightly closed. Store in a well-ventilated place. Keep cool. Store locked up. Dispose of contents/container in accordance with local/regional/national/international regulations.

Classification system:

NFPA ratings (scale 0 - 4)



HMIS-ratings (scale 0 - 4)



Other hazards

Results of PBT and vPvB assessment

- PBT: Not applicable. vPvB: Not applicable.

3 Composition/information on ingredients

Chemical characterization: Mixtures

Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: CAS number and component name, and 2 rows: CAS: 64-17-5 ethanol (≤25%), CAS: 67-63-0 propan-2-ol (≤25%)

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 2)

4 First-aid measures

- Description of first aid measures
After inhalation: Supply fresh air; consult doctor in case of complaints.
After skin contact: Generally the product does not irritate the skin.
After eye contact: Rinse opened eye for several minutes under running water.
After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed
Indication of any immediate medical attention and special treatment needed

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: CO2, extinguishing powder or water spray.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures
Environmental precautions: Dilute with plenty of water.
Methods and material for containment and cleaning up: Absorb with liquid-binding material.
Reference to other sections
Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, chemical name, and concentration (ppm). Rows include ethanol (1,800 ppm) and propan-2-ol (400 ppm).

PAC-2:

Table with 3 columns: CAS number, chemical name, and concentration (ppm). Rows include ethanol (3300* ppm) and propan-2-ol (2000* ppm).

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 3)

PAC-3:

Table with 3 columns: CAS number, chemical name, and limit value. Rows include ethanol (15000* ppm) and propan-2-ol (12000** ppm).

7 Handling and storage

- Handling:
Precautions for safe handling: Ensure good ventilation/exhaustion at the workplace.
Information about protection against explosions and fires: Keep ignition sources away - Do not smoke.
Conditions for safe storage, including any incompatibilities:
Storage:
Requirements to be met by storerooms and receptacles: Store in a cool location.
Information about storage in one common storage facility: Not required.
Further information about storage conditions: Keep receptacle tightly sealed.
Specific end use(s): Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters

Components with limit values that require monitoring at the workplace:

Table with 2 columns: CAS number and exposure limits. Rows include ethanol (CAS: 64-17-5) and propan-2-ol (CAS: 67-63-0) with PEL, REL, TLV, and BEI values.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 4)

Ingredients with biological limit values:

CAS: 67-63-0 propan-2-ol

Table with 2 columns: BEI, Value. BEI: 40 mg/L. Value: Medium: urine, Time: end of shift at end of workweek, Parameter: Acetone (background, nonspecific)

Additional information: The lists that were valid during the creation were used as basis.

Exposure controls

Personal protective equipment:

General protective and hygienic measures:

- Keep away from foodstuffs, beverages and feed. Immediately remove all soiled and contaminated clothing. Wash hands before breaks and at the end of work. Avoid contact with the eyes. Avoid contact with the eyes and skin.

Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection:



Tightly sealed goggles

US

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 5)

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

Form: Liquid
Color: Clear
Odor: Alcohol-like
Odor threshold: Not determined.

pH-value: Not determined.

Change in condition

Melting point/Melting range: Undetermined.
Boiling point/Boiling range: 78 °C (172 °F)

Flash point: 13 °C (55 °F)

Flammability (solid, gaseous): Not applicable.

Ignition temperature: 425 °C (797 °F)

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product is not explosive. However, formation of explosive air/vapor mixtures are possible.

Explosion limits:

Lower: 2.0 Vol %
Upper: 15.0 Vol %

Vapor pressure at 20 °C (68 °F): 59 hPa (44 mm Hg)

Density: Not determined.
Relative density: Not determined.
Vapor density: Not determined.
Evaporation rate: Not determined.

Solubility in / Miscibility with

Water: Fully miscible.

Partition coefficient (n-octanol/water): Not determined.

Viscosity:

Dynamic: Not determined.
Kinematic: Not determined.

Solvent content:

Organic solvents: 50.0 %
VOC content: 50.0 %
500.0 g/l / 4.17 lb/gl

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 6)

Other information No further relevant information available.

10 Stability and reactivity

- Reactivity No further relevant information available.
Chemical stability This product is normally stable under anticipated circumstances of use and storage.
Thermal decomposition / conditions to be avoided: Products of thermal decomposition of this material would include carbon monoxide and carbon dioxide.
Possibility of hazardous reactions No dangerous reactions known.
Conditions to avoid No further relevant information available.
Incompatible materials: Acids and strong oxidizers
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects May be harmful by inhalation. Material is irritating to mucous membranes and upper respiratory tract.
Acute toxicity:
Primary irritant effect:
on the skin: No irritant effect.
on the eye: Causes severe eye irritation Irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information: Inhalation of vapors/mists, and sprays of this product can cause dizziness, drowsiness and other central nervous system effects. The product shows the following dangers according to internally approved calculation methods for preparations: Irritant

Carcinogenic categories

IARC (International Agency for Research on Cancer)

CAS: 67-63-0 propan-2-ol 3

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 7)

- Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes:
Water hazard class 1 (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation:
Dispose of contents in accordance with local/regional/national, and international recommendations.
Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Label and Value. Rows include UN-Number (UN1993), UN proper shipping name (Flammable liquids, n.o.s. (Isopropanol, Ethanol)), Transport hazard class(es) (3), and Class (3 Flammable liquids). Includes a Flammable Liquid hazard pictogram.

(Contd. on page 9)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 8)

Table with hazard information including Label (3), IMDG/IATA (Flammable liquids 3), Packing group (II), Environmental hazards (Not applicable), Special precautions (Warning: Flammable liquids), and UN Model Regulation (UN 1993 FLAMMABLE LIQUIDS, N.O.S. (ISOPROPANOL, ETHANOL), 3, II).

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture Sara

Section 355 (extremely hazardous substances):

None of the ingredients is listed.

Section 313 (Specific toxic chemical listings):

CAS: 67-63-0 propan-2-ol

TSCA (Toxic Substances Control Act):

All ingredients are listed.

(Contd. on page 10)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 9)

Proposition 65

Chemicals known to cause cancer:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

CAS: 64-17-5 ethanol A3

CAS: 67-63-0 propan-2-ol A4

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms GHS02, GHS07

Signal word Danger

Hazard-determining components of labeling:

propan-2-ol
ethanol

Hazard statements

Highly flammable liquid and vapor.
Causes serious eye irritation.
May cause drowsiness or dizziness.

Precautionary statements

Keep away from heat/sparks/open flames/hot surfaces. No smoking.
Use explosion-proof electrical/ventilating/lighting/equipment.
Avoid breathing dust/fume/gas/mist/vapors/spray
Wear protective gloves / eye protection / face protection.
Ground/bond container and receiving equipment.
Use only non-sparking tools.
Take precautionary measures against static discharge.
Wash thoroughly after handling.
Use only outdoors or in a well-ventilated area.
If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
Continue rinsing.
IF INHALED: Remove person to fresh air and keep comfortable for breathing.
Call a POISON CENTER/doctor if you feel unwell.
If eye irritation persists: Get medical advice/attention.

(Contd. on page 11)

**Safety Data Sheet**
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 11/10/2016

Trade name: ZymoBIOMICS MagWash 2

(Contd. of page 10)

In case of fire: Use for extinction: CO₂, powder or water spray.
Store in a well-ventilated place. Keep container tightly closed.
Store in a well-ventilated place. Keep cool.
Store locked up.
Dispose of contents/container in accordance with local/regional/national/international regulations.
· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com· **Date of preparation / last revision** 08/22/2017 / -· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
DOT: US Department of Transportation
IATA: International Air Transport Association
ACGIH: American Conference of Governmental Industrial Hygienists
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit
BEI: Biological Exposure Limit
Flam. Liq. 2: Flammable liquids – Category 2
Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A
STOT SE 3: Specific target organ toxicity (single exposure) – Category 3



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 06/20/2016

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS DNase/RNase Free Water
- **Article number:** D4302-5-1, D4302-5-10, D4032-5-50
- **CAS Number:**
7732-18-5
- **EC number:**
231-791-2
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The substance is not classified according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

· **vPvB:** Not applicable.

(Contd. of page 1)

3 Composition/information on ingredients

- **Chemical characterization: Substances**
- **CAS No. Description**
7732-18-5 water, distilled, conductivity or of similar purity
- **Identification number(s)**
- **EC number:** 231-791-2

4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· **PAC-1:**
Substance is not listed.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 2)

· PAC-2:

Substance is not listed.

· PAC-3:

Substance is not listed.

7 Handling and storage

- **Handling:**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:** Not required.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Not required.
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Goggles recommended during refilling.



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 3)

9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· Form:	Liquid
· Color:	Clear
· Odor:	Odorless
· Odor threshold:	Not determined.

· **pH-value:** Not determined.

· **Change in condition**

· Melting point/Melting range:	0 °C (32 °F)
· Boiling point/Boiling range:	100 °C (212 °F)

· **Flash point:** Not applicable.

· **Flammability (solid, gaseous):** Not applicable.

· **Ignition temperature:**

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Not determined.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

· Lower:	Not Applicable
· Upper:	Not Applicable

· **Vapor pressure at 20 °C (68 °F):** 23 hPa (17 mm Hg)

· **Density at 20 °C (68 °F):** 1 g/cm³ (8.345 lbs/gal)

· **Relative density** Not determined.

· **Vapor density** Not determined.

· **Evaporation rate** Not determined.

· **Solubility in / Miscibility with**

· **Water:** Fully miscible.

· **Partition coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

· Dynamic at 20 °C (68 °F):	0.952 mPas
· Kinematic:	Not determined.

· **Solvent content:**

· Organic solvents:	0.0 %
· Water:	100.0 %

· **VOC content:** 0.0 g/l / 0.00 lb/gl

· **Other information** No further relevant information available.



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 4)

10 Stability and reactivity

- Reactivity No further relevant information available.
Chemical stability
Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
Possibility of hazardous reactions No dangerous reactions known.
Conditions to avoid No further relevant information available.
Incompatible materials: No further relevant information available.
Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
Acute toxicity:
Primary irritant effect:
on the skin: No irritant effect.
on the eye: No irritating effect.
Sensitization: No sensitizing effects known.
Additional toxicological information:
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.
The substance is not subject to classification.

Carcinogenic categories

IARC (International Agency for Research on Cancer)

Substance is not listed.

NTP (National Toxicology Program)

Substance is not listed.

OSHA-Ca (Occupational Safety & Health Administration)

Substance is not listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes: Generally not hazardous for water
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

(Contd. on page 6)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 5)

Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Smaller quantities can be disposed of with household waste.
Uncleaned packagings:
Recommendation: Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

Table with 2 columns: Regulatory information and Status. Rows include UN-Number, proper shipping name, transport hazard class(es), packing group, environmental hazards, special precautions for user, transport in bulk, and model regulation.

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Table with 2 columns: Regulatory section and Status. Rows include Section 355 (extremely hazardous substances), Section 313 (Specific toxic chemical listings), and TSCA (Toxic Substances Control Act).

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 6)

Proposition 65

Chemicals known to cause cancer:

Substance is not listed.

Chemicals known to cause reproductive toxicity for females:

Substance is not listed.

Chemicals known to cause reproductive toxicity for males:

Substance is not listed.

Chemicals known to cause developmental toxicity:

Substance is not listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

Substance is not listed.

TLV (Threshold Limit Value established by ACGIH)

Substance is not listed.

NIOSH-Ca (National Institute for Occupational Safety and Health)

Substance is not listed.

GHS label elements Void

Hazard pictograms Void

Signal word Void

Hazard statements Void

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Department issuing SDS:

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

Contact: sds@zymoresearch.com

Date of preparation / last revision 08/22/2017 / -

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 06/20/2016

Trade name: ZymoBIOMICS DNase/RNase Free Water

(Contd. of page 7)

HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

1 Identification

- **Product identifier**
- **Trade name:** ZymoBIOMICS™ MagBinding Beads
- **Article number:** D4302-6-12
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product safety department
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS07

- Acute Tox. 4 H302 Harmful if swallowed.
- Skin Irrit. 2 H315 Causes skin irritation.
- Eye Irrit. 2A H319 Causes serious eye irritation.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Hazard-determining components of labeling:**
guanidinium chloride
- **Hazard statements**
Harmful if swallowed.
Causes skin irritation.
Causes serious eye irritation.
- **Precautionary statements**
Wear protective gloves / eye protection / face protection.
Wash thoroughly after handling.
Do not eat, drink or smoke when using this product.
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
Continue rinsing.
Specific treatment (see on this label).
IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
If skin irritation occurs: Get medical advice/attention.

(Contd. on page 2)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 1)

If eye irritation persists: Get medical advice/attention.
Rinse mouth.
IF ON SKIN: Wash with plenty of water.
Take off contaminated clothing and wash it before reuse.
Dispose of contents/container in accordance with local/regional/national/international regulations.

- Classification system:
NFPA ratings (scale 0 - 4)



- HMIS-ratings (scale 0 - 4)

Table with 3 rows: HEALTH 1, FIRE 0, REACTIVITY 0

- Other hazards
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

3 Composition/information on ingredients

- Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: CAS: 50-01-1 guanidinium chloride, ≤100%

4 First-aid measures

- Description of first aid measures
General information: Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.
After inhalation: In case of unconsciousness place patient stably in side position for transportation.
After skin contact: Immediately wash with water and soap and rinse thoroughly.
After eye contact: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
After swallowing: Do not induce vomiting; immediately call for medical help.
Information for doctor:
Most important symptoms and effects, both acute and delayed No further relevant information available.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 2)

- Indication of any immediate medical attention and special treatment needed
No further relevant information available.

5 Fire-fighting measures

- Extinguishing media
Suitable extinguishing agents: CO2, extinguishing powder or water spray.
Advice for firefighters
Protective equipment: Wear protective clothing.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures
Environmental precautions: Do not allow to enter sewers/ surface or ground water.
Methods and material for containment and cleaning up:
Reference to other sections
Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, chemical name, and concentration (mg/m3). Rows include guanidinium chloride, diiron trioxide, and silicon dioxide.

PAC-2:

Table with 3 columns: CAS number, chemical name, and concentration (mg/m3). Rows include guanidinium chloride, diiron trioxide, and silicon dioxide.

PAC-3:

Table with 3 columns: CAS number, chemical name, and concentration (mg/m3). Rows include guanidinium chloride, diiron trioxide, and silicon dioxide.

7 Handling and storage

- Handling:
Precautions for safe handling
Ensure good ventilation/exhaustion at the workplace.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 3)

- Prevent formation of aerosols.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.
- **Breathing equipment:**
In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.
- **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 4)

· Eye protection:



Tightly sealed goggles

9 Physical and chemical properties

· Information on basic physical and chemical properties

· General Information

· Appearance:

- Form: Suspension
Color: Dark
Odor: No information available
Odor threshold: Not determined.

· pH-value: Not determined.

· Change in condition

- Melting point/Melting range: Undetermined.
Boiling point/Boiling range: Undetermined.

· Flash point: Not applicable.

· Flammability (solid, gaseous): Not applicable.

· Ignition temperature: > 370 °C (> 698 °F)

· Decomposition temperature: Not determined.

· Auto igniting: Product is not selfigniting.

· Danger of explosion: Product does not present an explosion hazard.

· Explosion limits:

- Lower: Not Applicable
Upper: Not Applicable

· Vapor pressure: Not determined.

- Density: Not determined.
Relative density: Not determined.
Vapor density: Not determined.
Evaporation rate: Not determined.

· Solubility in / Miscibility with

Water: Not miscible or difficult to mix.

· Partition coefficient (n-octanol/water): Not determined.

· Viscosity:

Dynamic: Not determined.

(Contd. on page 6)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 5)

Table with 2 columns: Property (Kinematic, Solvent content, Organic solvents, VOC content, Solids content, Other information) and Value (Not determined, 0.0 %, 0.0 g/l / 0.00 lb/gl, 80.0 %, No further relevant information available).

10 Stability and reactivity

- Reactivity No further relevant information available.
· Chemical stability
· Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
· Possibility of hazardous reactions No dangerous reactions known.
· Conditions to avoid No further relevant information available.
· Incompatible materials: No further relevant information available.
· Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on toxicological effects
· Acute toxicity:

Table with 2 columns: LD/LC50 values that are relevant for classification and CAS: 50-01-1 guanidinium chloride. Row 1: Oral LD50 475 mg/kg (rat)

- Primary irritant effect:
· on the skin: Irritant to skin and mucous membranes.
· on the eye: Irritating effect.
· Sensitization: No sensitizing effects known.

Additional toxicological information: The product shows the following dangers according to internally approved calculation methods for preparations: Harmful, Irritant

- Carcinogenic categories

Table with 3 columns: IARC (International Agency for Research on Cancer), CAS number, and hazard level. Row 1: CAS: 1309-37-1 diiron trioxide, 3. Row 2: CAS: 7631-86-9 silicon dioxide, chemically prepared, 3.

- NTP (National Toxicology Program)

None of the ingredients is listed.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 6)

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
· Aquatic toxicity: No further relevant information available.
· Persistence and degradability No further relevant information available.
· Behavior in environmental systems:
· Bioaccumulative potential No further relevant information available.
· Mobility in soil No further relevant information available.
· Additional ecological information:
· General notes:
Water hazard class 1 (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
· Results of PBT and vPvB assessment
· PBT: Not applicable.
· vPvB: Not applicable.
· Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
· Recommendation:
Dispose of contents in accordance with local/regional/national, and international recommendations.
· Uncleaned packagings:
· Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.

14 Transport information

- UN-Number
· DOT, ADN, IMDG, IATA not regulated
· UN proper shipping name
· DOT, ADN, IMDG, IATA not regulated
· Transport hazard class(es)
· DOT, ADN, IMDG, IATA
· Class not regulated
· Packing group
· DOT, IMDG, IATA not regulated

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 7)

- Environmental hazards: Not applicable.
Special precautions for user: Not applicable.
Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: Not applicable.
UN Model Regulation: not regulated

15 Regulatory information

- Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Section 355 (extremely hazardous substances):

None of the ingredients is listed.

Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

TSCA (Toxic Substances Control Act):

All ingredients are listed.

Proposition 65

Chemicals known to cause cancer:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

CAS: 1309-37-1 diiron trioxide

A4

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements The product is classified and labeled according to the Globally Harmonized System (GHS).

Hazard pictograms GHS07

Signal word Warning

Hazard-determining components of labeling:

guanidinium chloride

(Contd. on page 9)



Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 8)

· **Hazard statements**

Harmful if swallowed.
Causes skin irritation.
Causes serious eye irritation.

· **Precautionary statements**

Wear protective gloves / eye protection / face protection.
Wash thoroughly after handling.
Do not eat, drink or smoke when using this product.
If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
Continue rinsing.
Specific treatment (see on this label).
IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
If skin irritation occurs: Get medical advice/attention.
If eye irritation persists: Get medical advice/attention.
Rinse mouth.
IF ON SKIN: Wash with plenty of water.
Take off contaminated clothing and wash it before reuse.
Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 08/22/2017 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
DOT: US Department of Transportation
IATA: International Air Transport Association
ACGIH: American Conference of Governmental Industrial Hygienists
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
LC50: Lethal concentration, 50 percent
LD50: Lethal dose, 50 percent
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety

(Contd. on page 10)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 10/10

Safety Data Sheet
acc. to OSHA HCS

Printing date 08/22/2017

Reviewed on 03/24/2017

Trade name: ZymoBIOMICS™ MagBinding Beads

(Contd. of page 9)

OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit
Acute Tox. 4: Acute toxicity – Category 4
Skin Irrit. 2: Skin corrosion/irritation – Category 2
Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A

US