



# ZR DNA Sequencing Clean-up Kits™

- **Long Read Lengths:** Complete elimination of “dye blobs” for high-quality Phred scores.
- **Reusable:** Columns can be reused.
- **Flexible Elutions:** 6-20 µl elution volumes allow for direct loading of samples without precipitation or drying.

## Simple, Innovative Workflow



## Ultra-pure DNA without “Dye Blobs”



Sequencing chromatogram of pGEM® DNA generated using an ABI 3730xl DNA analyzer. DNA was labeled with ABI BigDye® v3.1 Terminators and cleaned using the ZR DNA Sequencing Clean-up Kit™.

| Product                             | Cat. No. | Size         |
|-------------------------------------|----------|--------------|
| ZR DNA Sequencing Clean-Up Kits™    | D4050    | 50 preps     |
|                                     | D4051    | 200 preps    |
| ZR-96 DNA Sequencing Clean-Up Kits™ | D4052    | 2 x 96 preps |
|                                     | D4053    | 4 x 96 preps |

