

Safety Data Sheet acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

1 Identification

· **Product identifier**

· **Trade name:** DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

· **Article number:**

R1100-50, R1100-250, R1100-8-S, R1200-25, R1200-125, R1150, R1101, R1102, R1103, R1104, R1105, R1106, R1107, R1108, R1109, R1210

· **Application of the substance / the mixture** Laboratory Reagent

· **Details of the supplier of the safety data sheet**

· **Manufacturer/Supplier:**

Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com

· **Information department:** Product safety department

· **Emergency telephone number:**

During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

· **Classification of the substance or mixture**



GHS07

Acute Tox. 4 H302 Harmful if swallowed.

Skin Irrit. 2 H315 Causes skin irritation.

Eye Irrit. 2B H320 Causes eye irritation.

· **Label elements**

· **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms** GHS07

· **Signal word** Warning

· **Hazard statements**

Harmful if swallowed.
Causes skin and eye irritation.

· **Precautionary statements**

Wash thoroughly after handling.
Do not eat, drink or smoke when using this product.
Wear protective gloves.
If swallowed: Call a poison center/doctor if you feel unwell.

(Contd. on page 2)



Safety Data Sheet

acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 1)

If on skin: Wash with plenty of water.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Specific treatment (see on this label).

Rinse mouth.

If skin irritation occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

Take off contaminated clothing and wash it before reuse.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



Health = 2

Fire = 0

Reactivity = 0

· **HMIS-ratings (scale 0 - 4)**



Health = 2

Fire = 0

Reactivity = 0

· **Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** This product is a proprietary solution.

· **Dangerous components:** Void

4 First-aid measures

· **Description of first aid measures**

· **After inhalation:** Supply fresh air; consult doctor in case of complaints.

· **After skin contact:** Immediately wash with water and soap and rinse thoroughly.

· **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.

· **After swallowing:** Do not induce vomiting; immediately call for medical help.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

(Contd. on page 3)



Safety Data Sheet

acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

Trade name: *DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit*

(Contd. of page 2)

· Indication of any immediate medical attention and special treatment needed

No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
Not considered a fire hazard. Use extinguishing media appropriate surrounding fire conditions.
Use fire fighting measures that suit the environment.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**
Ensure adequate ventilation.
Dispose contaminated material as waste according to item 13.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and storage

- **Handling:**
- **Precautions for safe handling** No special precautions are necessary if used correctly.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s)** Laboratory reagent

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.

(Contd. on page 4)



Safety Data Sheet

acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 3)

- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:** Wash hands before breaks and at the end of work.
- **Breathing equipment:** Not required.
- **Protection of hands:** Use appropriate gloves. Avoid contact with skin.
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Avoid contact with eyes.
- **Body protection:**
Avoid contact with eyes, skin, and clothing. Wash thoroughly after handling. Remove contaminated clothing and wash before use. Avoid ingestion and inhalation.

9 Physical and chemical properties

- **Information on basic physical and chemical properties**
- **General Information**
- **Appearance:**

Form:	Liquid
Color:	Colorless
Odor:	No information available
Odor threshold:	Not determined.
- **pH-value:** 5.0-7.0
- **Change in condition**

Melting point/Melting range:	Undetermined.
Boiling point/Boiling range:	Undetermined.
- **Flash point:** Not applicable.
- **Flammability (solid, gaseous):** Not applicable.
- **Ignition temperature:**
- **Decomposition temperature:** Not determined.
- **Auto igniting:** Product is not selfigniting.
- **Danger of explosion:** Product does not present an explosion hazard.

(Contd. on page 5)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 5/9

Safety Data Sheet

acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 4)

· Explosion limits:	
Lower:	Not Applicable
Upper:	Not Applicable
· Vapor pressure:	Not determined.
· Density:	Not determined.
· Relative density	Not determined.
· Vapor density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with Water:	Fully miscible.
· Partition coefficient (n-octanol/water):	Not determined.
· Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined.
· Solvent content:	
VOC content:	0,00 %
Solids content:	0,0 %
· Other information	No further relevant information available.

10 Stability and reactivity

- **Reactivity** Strong oxidizing agents and acids.
- **Possibility of hazardous reactions** Strong oxidizing agents and acids.
- **Conditions to avoid** Strong oxidizing agents and acids.
- **Incompatible materials:** Strong oxidizing agents and acids.

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the eye:** Irritating to eyes.
- **Additional toxicological information:**
- **Carcinogenic categories**
- **IARC (International Agency for Research on Cancer)**

None of the ingredients is listed.

(Contd. on page 6)

US



Safety Data Sheet

acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 5)

· **NTP (National Toxicology Program)**

None of the ingredients is listed.

· **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**
Dispose of contents in accordance with local/regional/national, and international recommendations.
- **Uncleaned packagings:**
- **Recommendation:**
Dispose of container in accordance with local/regional/national and international recommendations.

14 Transport information

- | | |
|-------------------------------------|----------------------|
| · UN-Number | Not dangerous goods. |
| · DOT, ADN, IMDG, IATA | not regulated |
| · UN proper shipping name | Not dangerous goods. |
| · DOT, ADN, IMDG, IATA | not regulated |
| · Transport hazard class(es) | Not dangerous goods. |
| · DOT, ADN, IMDG, IATA | |
| · Class | not regulated |
| · Packing group | Not dangerous goods. |

(Contd. on page 7)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 7/9

Safety Data Sheet

acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

Trade name: *DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit*

(Contd. of page 6)

· DOT, IMDG, IATA	not regulated
· Environmental hazards:	Not applicable.
· Special precautions for user	Not applicable.
· Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code	Not applicable.
· UN "Model Regulation":	not regulated

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

- **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

- **Section 313 (Specific toxic chemical listings):**

None of the ingredients is listed.

- **TSCA (Toxic Substances Control Act):**

None of the ingredients is listed.

- **Proposition 65**

- **Chemicals known to cause cancer:**

None of the ingredients is listed.

- **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

- **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

- **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

- **Carcinogenic categories**

- **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

- **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

- **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

(Contd. on page 8)



Safety Data Sheet

acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

Trade name: *DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit*

(Contd. of page 7)

- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Hazard statements**
 - Harmful if swallowed.
 - Causes skin and eye irritation.
- **Precautionary statements**
 - Wash thoroughly after handling.
 - Do not eat, drink or smoke when using this product.
 - Wear protective gloves.
 - If swallowed: Call a poison center/doctor if you feel unwell.
 - If on skin: Wash with plenty of water.
 - If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 - Specific treatment (see on this label).
 - Rinse mouth.
 - If skin irritation occurs: Get medical advice/attention.
 - If eye irritation persists: Get medical advice/attention.
 - Take off contaminated clothing and wash it before reuse.
 - Dispose of contents/container in accordance with local/regional/national/international regulations.
- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Department issuing SDS:**
 - Zymo Research Corp.
 - Safety Department
 - 17062 Murphy Ave.
 - Irvine, CA 92614
 - USA
 - Phone: 1-949-679-1190 or 1-888-882-9682
- **Contact:** sds@zymoresearch.com
- **Date of preparation / last revision** 04/26/2018 / -
- **Abbreviations and acronyms:**
 - ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
 - IMDG: International Maritime Code for Dangerous Goods
 - DOT: US Department of Transportation
 - IATA: International Air Transport Association
 - GHS: Globally Harmonised System of Classification and Labelling of Chemicals
 - EINECS: European Inventory of Existing Commercial Chemical Substances
 - ELINCS: European List of Notified Chemical Substances
 - CAS: Chemical Abstracts Service (division of the American Chemical Society)
 - VOC: Volatile Organic Compounds (USA, EU)
 - PBT: Persistent, Bioaccumulative and Toxic

(Contd. on page 9)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 9/9

Safety Data Sheet
acc. to OSHA HCS

Printing date 04/26/2018

Reviewed on 05/06/2016

Trade name: DNA/RNA Shield - 50 ml, 250 ml & sample; DNA/RNA Shield (2X Concentrate) - 25 ml & 125 ml; DNA/RNA Shield Blood Collection Tube, DNA/RNA Shield Fecal Collection Tube, DNA/RNA Shield Collection Tube, DNA/RNA Shield Lysis Tube (Microbe), DNA/RNA Shield Lysis Tube w/ Swab (Microbe), DNA/RNA Shield Lysis Tube (Tissue); DNA/RNA Shield Collection Tube w/ Swab; DNA/RNA Shield Saliva Collection Kit

(Contd. of page 8)

vPvB: very Persistent and very Bioaccumulative
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit
Acute Tox. 4: Acute toxicity – Category 4
Skin Irrit. 2: Skin corrosion/irritation – Category 2
Eye Irrit. 2B: Serious eye damage/eye irritation – Category 2B

US