

www.biolinkk.com





Cellulose Acetate Membrane Electrophoresis (BL-HE7)

The Biolinkk Cellulose Acetate Membrane Electrophoresis System (Model BL-HE7) features a double bridge design for enhanced performance, a cooling device for temperature control, and adjustable rods for flexible electrophoresis distances. Ideal for precise and efficient electrophoresis applications.



Highlights

- Double bridge design for enhanced flexibility and performance.
- Integrated cooling device for constant temperature control.
- Adjustable rods to customize electrophoretic distance based on medium length.

The Biolinkk Cellulose Acetate Membrane Electrophoresis System (Model BL-HE7) offers advanced capabilities for high-precision electrophoresis. Featuring a double bridge style with bridge sizes of 70x250mm and 90x250mm, the system accommodates different membrane sizes (70x90mm) and maintains consistent results with a large buffer volume of 1000 ml. Equipped with a cooling device, it ensures constant temperature throughout the electrophoresis process, making it ideal for temperature-sensitive experiments. The adjustable rod allows the user to customize the electrophoretic distance based on the medium length, providing flexibility for various applications. Maintenance is simplified with unfolded electrodes for easy repair and part replacement. For added safety, the system includes an auto-switch-off function when the lid is removed, ensuring secure operation.



Features

- Bridge Size: 70x250mmx2; 90x250mmx2.
- Membrane Size: 70x90mm.
- Buffer Volume: 1000 ml.
- Double Bridge Design: Enhances flexibility and efficiency during electrophoresis.
- Cooling Device: Maintains constant temperature for sensitive samples.
- Adjustable Rods: Customize electrophoresis distance to meet experimental requirements.
- Unfolded Electrodes: Simplified maintenance with easy part replacement.
- Safety Feature: Auto-switch-off when the lid is removed for secure operation.



| Model | BL-HE7 |
|-------------------|-------------------------|
| Bridge Size (WxL) | 70x250mmx2; 90x250mmx2. |
| Buffer Volume | 1000ml |
| Dimension (LxWxH) | 350×260×65 (mm) |
| Net weight | 1.6 (kg) |



Add: Third Floor, Plot No. 25, Community Centre, Block B-1, Janakpuri, Delhi-110058

To order contact us:

Email: customerservice@biolinkk.com Tel: +91-11-20003-5753 / 25553658 www.bolinkk.com













