




Microbiology Tools for Research




BioVision is excited to introduce numerous staining kits and reagents for pathology including the bacterial staining kits namely Acid-Fast bacteria stain kit, Fite's stain kit, Gram stain kit; Mucicarmine staining kit for staining mucins; Wright-Giemsa stain kit for differential staining of blood smears and many more staining kits and reagents. These special staining kits are intended for use in the histological visualization of organisms, cell types, structures in tissue sections, cell smears etc. These kits play a critical role in tissue or cell-based diagnosis or research. Additionally, BioVision offers many routine counterstains, kits and reagents including hematoxylin, eosin, modified silver stain etc.

Applications:

- **Histological Applications**
- **For *in vitro* Diagnostic Use only**

Key Features:

- **Rapid procedure**
- **Many staining options to choose from**
- **Can be used on fixed paraffin embedded, formalin fixed, frozen tissue sections or cell smears**



Staining using **A.** Acid Fast Bacteria stain kit. **B.** Fite's stain kit. **C.** Mucicarmine stain kit. **D.** Wright-Giemsa stain kit

155 S. Milipitas Blvd, Milipitas, CA 95035

T: 408-493-1800 F: 408-493-1801

Toll Free: 800-891-9699 (US Only)

BioVision
BioVision Incorporated



Check out the full list of newly added microbiology staining kits and reagents available from BioVision

Product Name	Cat. No.	Size
Acid Fast Bacteria Stain Kit	K1424	30 ml, 125 ml
Fite's Stain Kit	K1425	125 ml
Geimsa Stain Kit	K1426	30 ml, 500 ml
GMS Stain Kit	K1427	125 ml
Gram Stain Kit	K1428	30 ml, 125 ml
Modified Gram Stain Kit	K1429	30 ml, 125 ml
Modified Southgate's Mucicarmine Stain Kit	K1430	30 ml, 125 ml
Orcein Stain Kit	K1432	30 ml, 125 ml
Papanicolaou Stain Kit	K1440	30 ml, 500 ml
Periodic Acid Schiff Diastase Stain Kit	K1433	30 ml, 250 ml
Pneumocystis Stain Kit	K1434	30 ml, 125 ml
Southgate's Mucicarmine Solution	M1430	125 ml
Steiner Stain Kit	K1435	125 ml
Twort's Counterstain Kit	K1436	125 ml
Warthin-Starry Stain Kit	K1437	125 ml
Wright-Giemsa Solution	M1438	500 ml, 1000 ml
Wright-Giemsa Stain Kit	K1438	30 ml, 500 ml

Please visit <https://www.biovision.com/products/microbiology/staining-kits.html> to see the full range of Microbiology Staining Kits from BioVision and our Microbiology Staining Kits Flyer.