


# Monoclonal Secondary Antibodies

SPECIAL OFFER!

biolinkk  
Innovation to support science...



Terms :

Offer is Valid up to 31<sup>st</sup> October 2019

| Cat No.                             | Description                                 | Pack Size | List Price | Offer Price | Cat No.                           | Description                   | Pack Size | List Price | Offer Price |
|-------------------------------------|---|-----------|------------|-------------|-----------------------------------|-------------------------------|-----------|------------|-------------|
| <b>HRP Conjugated</b>               |   |           |            |             | <b>Biotin Conjugated</b>          |                               |           |            |             |
| 7100-05                             | Streptavidin                                | 1.0mL     | 15,840.00  | 14,256.00   | 1036-08                           | Goat Anti-Mouse IgG(H+L)      | 1.0mg     | 10,368.00  | 9,331.00    |
| 1033-05                             | Goat Anti-Mouse IgG Fc                      | 1.0mL     | 14,400.00  | 12,960.00   | 4020-08                           | Goat Anti-Rabbit IgM          | 1.0mg     | 14,400.00  | 12,960.00   |
| 1030-05                             | Goat Anti-Mouse IgG, Human ads              | 1.0mL     | 15,840.00  | 14,256.00   | 6420-08                           | Donkey Anti-Goat IgG(H+L)     | 1.0mg     | 12,360.00  | 11,340.00   |
| 4030-05                             | Goat Anti-Rabbit IgG                        | 1.0mL     | 11,520.00  | 10,886.00   | 6022-08                           | Rabbit Anti-Goat IgG          | 0.5mg     | 21,600.00  | 19,440.00   |
| 6420-05                             | Donkey Anti-Goat IgG(H+L), Rat/Mouse SP ads | 1.0mL     | 15,840.00  | 14,256.00   | 2048-08                           | Goat Anti Human IGG           | 0.5mg     | 11,232.00  | 10,108.00   |
| 6160-05                             | Rabbit Anti-Goat IgG(H+L)                   | 1.0mL     | 12,614.00  | 11,818.00   | 3050-08                           | Goat Anti-Rat IgG(H+L)        | 1.0mg     | 19,800.00  | 17,820.00   |
| 6050-05                             | Goat Anti-Porcine IgG(H+L)                  | 1.0mL     | 12,096.00  | 10,886.00   | 6411-08                           | Donkey Anti-Mouse IgG(H+L)    | 1.0mg     | 10,200.00  | 9180.00     |
| 6100-05                             | Goat Anti-Chicken IgY(H+L)                  | 1.0mL     | 12,086.00  | 10,886.00   | 6441-08                           | Donkey Anti-Rabbit IgG(H+L)   | 1.0mg     | 7,200.00   | 6,480.00    |
| 4060-05                             | Mouse Anti-Rabbit Light Chain               | 1.0mL     | 39,600.00  | 35,640.00   | <b>Texas Red Conjugated</b>       |                               |           |            |             |
| 6411-05                             | Donkey Anti-Mouse IgG(H+L)                  | 1.0mL     | 11,520.00  | 10,368.00   | 1030-07                           | Goat Anti-Mouse IgG           | 1.0mg     | 12,600.00  | 11,340.00   |
| 6441-05                             | Donkey Anti-Rabbit IgG(H+L)                 | 1.0mL     | 11,520.00  | 10,368.00   | 4050-07                           | Goat Anti-Rabbit IgG(H+L)     | 1.0mg     | 10,224.00  | 9590.00     |
| 6090-05                             | Goat Anti-Guinea Pig IgG(H+L)               | 1.0mL     | 12,096.00  | 10,886.00   | 2010-07                           | Goat Anti-Human Ig            | 1.0mg     | 8,280.00   | 7,560.00    |
| 6040-05                             | Goat Anti-Equine IgG(H+L)                   | 1.0mL     | 12,096.00  | 10,886.00   | <b>Alexa Fluor 488 Conjugated</b> |                               |           |            |             |
| 4700-05                             | Mouse Anti-Monkey IgG                       | 1.0mL     | 28,800.00  | 25,920.00   | 1030-30                           | Goat Anti-Mouse IgG           | 1.0mg     | 18,000.00  | 16,200.00   |
| 3030-05                             | Goat Anti-Rat IgG                           | 1.0mL     | 14,112.00  | 13219.00    | 4030-30                           | Goat Anti-Rabbit IgG          | 1.0mg     | 15,000.00  | 13,500.00   |
| 6140-05                             | Rabbit Anti-Human IgG(H+L)                  | 1.0mL     | 11,232.00  | 10,108.00   | 6411-30                           | Donkey Anti-Mouse IgG(H+L)    | 1.0mg     | 18,600.00  | 16,740.00   |
| 6150-05                             | Rabbit Anti-Sheep IgG(H+L)                  | 1.0mL     | 12,614.00  | 11,818.00   | <b>Alexa Flour 647</b>            |                               |           |            |             |
| 2040-05                             | Goat Anti-Human IgG                         | 1.0mL     | 12,096.00  | 10,886.00   | 1030-31                           | Goat Anti-Mouse IgG           | 1.0mg     | 18,000.00  | 16,200.00   |
| 7300-05                             | Protein A                                   | 1.0mL     | 13,680.00  | 12,312.00   | 4010-31                           | Goat Anti-Rabbit Ig           | 1.0mg     | 16,200.00  | 14,580.00   |
| <b>Affinity-Purified Antibodies</b> |   |           |            |             | 6411-31                           | Donkey Anti-Mouse IgG(H+L)    | 1.0mg     | 18,600.00  | 16,740.00   |
| 1030-04                             | Goat Anti-Mouse IgG                         | 1.0mL     | 17,400.00  | 15,660.00   | <b>Rhodamine Conjugated</b>       |                               |           |            |             |
| 4020-04                             | Goat Anti-Rabbit IgM                        | 1.0mL     | 14,400.00  | 12,960.00   | 1030-03                           | Goat Anti-Mouse IgG           | 1.0mg     | 12,000.00  | 10,800.00   |
| 6420-04                             | Donkey Anti-Goat IgG(H+L), Rat/Mouse SP ads | 1.0mL     | 17,400.00  | 15,660.00   | 4030-03                           | Goat Anti-Rabbit IgG          | 1.0mg     | 10,944.00  | 10,238.00   |
| 6411-04                             | Donkey Anti-Mouse IgG(H+L)                  | 1.0mL     | 11,400.00  | 10,260.00   | 6160-03                           | Rabbit Anti-Goat IgG(H+L)     | 1.0mg     | 6,240.00   | 5,832.00    |
| 6441-04                             | Donkey Anti-Rabbit IgG(H+L)                 | 1.0mL     | 11,400.00  | 10,260.00   | 3020-03                           | Goat Anti-Rat IgM             | 1.0mg     | 10,800.00  | 9,720.00    |
| 6430-04                             | Donkey Anti-Rat IgG(H+L), Mouse SP ads      | 1.0mL     | 17,400.00  | 15,660.00   | 3030-03                           | Goat Anti-Rat IgG             | 1.0mg     | 9,120.00   | 8,532.00    |
| 2085-04                             | Goat Anti-Human Ig Fab                      | 1.0mL     | 18,480.00  | 17,388.00   | 2040-03                           | Goat Anti Human IGG           | 1.0mg     | 10,200.00  | 9,180.00    |
| <b>FITC Conjugated</b>              |   |           |            |             | <b>PE Conjugated</b>              |                               |           |            |             |
| 1030-02                             | Goat Anti-Mouse IgG, Human ads              | 1.0mg     | 14,400.00  | 12,960.00   | 1030-09                           | Goat Anti-Mouse IgG           | 0.5mg     | 22,800.00  | 20,520.00   |
| 4030-02                             | Goat Anti-Rabbit IgG                        | 1.0mg     | 11,520.00  | 10,368.00   | 4030-09                           | Goat Anti-Rabbit IgG          | 0.5mg     | 23,400.00  | 21,060.00   |
| 6420-02                             | Donkey Anti-Goat IgG(H+L), Rat/Mouse SP ads | 1.0mg     | 15,120.00  | 13,608.00   | 6157-09                           | Mouse Anti-Goat               | 0.1mg     | 24,000.00  | 21,600.00   |
| 6164-02                             | Rabbit Anti-Goat IgG(H+L), Human SP ads     | 1.0mg     | 11,520.00  | 10,368.00   | 3030-09                           | Goat Anti-Rat IgG             | 0.5mg     | 25,800.00  | 23,220.00   |
| 6441-02                             | Donkey Anti-Rabbit IgG(H+L)                 | 1.0mg     | 10,800.00  | 9,720.00    | 2040-09                           | Goat Anti-Human IgG           | 0.5mg     | 23,400.00  | 21,060.00   |
| 6411-02                             | Donkey Anti-Mouse IgG(H+L)                  | 1.0mg     | 10,800.00  | 9,720.00    | <b>Cy5 Conjugated</b>             |                               |           |            |             |
| 3030-02                             | Goat Anti-Rat IgG                           | 1.0mg     | 12,672.00  | 11,923.00   | 1030-15                           | Goat Anti-Mouse IgG           | 1.0mg     | 17,400.00  | 15,660.00   |
| 2040-02                             | Goat Anti-Human IgG                         | 1.0mg     | 12,240.00  | 11,016.00   | 4050-09                           | Cy5 Goat Anti-Rabbit IgG(H+L) | 0.5mg     | 38,400.00  | 34,560.00   |
|                                     |   |           |            |             | 3060-15                           | Mouse Anti-Rat IgG1           | 0.1mg     | 36,000.00  | 32,400.00   |