

YEAR END OFFER
UPTO 31st December

Special Offer







**THERMO SCIENTIFIC
 ART PIPETTE TIP
 are now Available**

- BARRIER



ART Pipette Tips

The ART Barrier provides an environmentally-responsible solution without compromising sterility. With offset wells that allow the tips to stack, this advanced system reduces waste and saves bench top and storage space.

Cat No.	Description	Pack	List Price	Offer Price
10 μ L				
2139	ART 10 MICROPOINT RACKED ST	96 Tips/Tray, 10 Trays/Pack	5,759.00	2,500.00
20 μ L				
2149	ART 20UL PIPET TIP RACKED ST	96 Tips/Tray, 10 Trays/Pack	5,670.00	2,800.00
200 μ L				
2069	ART 200 MICROPOINT RACKED ST	96 Tips/Tray, 10Trays/Pack	5,681.00	2,500.00
1000 μ L				
2279	ART 1000 SOFTFIT RACKED ST	96 Tips/Tray, 8 Trays/Pack	5,600.00	2,500.00