

Grain Hardness Meter



Model: BL-GHM-II

Ideal device for hardness of granular good such as rice, wheat etc...



SPECIAL OFFER

95,200/-

Characteristics :

1. Easy to operate, direct reading, high accuracy, elegant appearance.
2. Pressure head is made of stainless steel, which increases the service life.
3. Has two types: Pointer and Digital LED.

Technical Parameters

| | |
|------------------|--------|
| Test range | 20kg |
| Display Digital | LED |
| Maximum Diameter | 18mm |
| Stress | Manual |
| Accuracy | ±0.5% |
| Resolution | 0.1kg |

BIOLINKK

Building No. 1333 4th Floor, Office No: 412, Durga Chamber, Des Bandhu Gupta Road, Karol Bagh, New Delhi-110005
Web: www.biolinkk.com Tel: (+91)-011-45035753/011-28753658