

SEED TESTING INSTRUMENTS

Vacuum seed Counter & Placing apparatus

SPECIAL
OFFER
4,00,000



□ **Model : BL-VSCPA-50**

□ **Product introduction :**

Vacuum seed Counter is specially used for seed germination test in the process of several species, seed suction, seed placement. It mainly consists of three parts: workbench, extension pipette and suction head.

□ **Technical parameters**

- Sucker configuration: 5 pieces (1 large, 3 small, 1 round)
- Motor power: 1440W
- Working power supply: 220V 50Hz
- Power regulation range: 300W-1400W
- The aperture of 100 holes: 6mm/3mm/1.5mm/0.7mm

□ **Characteristics :**

- ◇ Adopt High-speed vacuum axial flow fan, with strong suction, and low noise.
- ◇ Built-in dirt isolation bag of the fuselage can easily remove the debris and moisture sucked into the fuselage.
- ◇ Extension pipette is connected to the plug-in handle at one end and has a vent through hole, which causes the seeds to fall down quickly and smoothly.
- ◇ Set the bottom of the bed table with pulleys and move flexibly. Equipped with a key control switch and 5 suction heads, which are of the same size and specification as the domestic common germination box. To meet the needs of different seed germination tests.

BIOLINKK

Building NO: 1333, 4th Floor, Office No:412, Durga Chamber, Desh Bandhu Gupta Road, Karol Bagh, New Delhi- 110005. (India),
Landline :- (+91)-11-45035753, (+91)-11-28753658, Website :- www.biolinkk.com