

DNA/RNA Shield Saliva Collection Kit

Cat no(s): R1210

All biological samples stabilized in DNA/RNA Shield™ are:

- (a) compatible with most workflows (e.g., Qiagen, Roche, Thermo-Fisher).
- (b) a direct input for high-throughput processing (e.g., Hamilton, Tecan, bioMérieux, PerkinElmer, Eppendorf, Promega, Opentrons).

Recommended Purification kits

	Cellular	Cell-Free
DNA	<p><i>Quick-DNA Plus Kits</i></p> <p>D4068 (Column) D4081 (Magbead)</p>	<p><i>Quick-cfDNA Serum & Plasma Kit</i></p> <p>D4076</p>
RNA	<p><i>Quick-RNA Miniprep Plus Kit</i></p> <p>R1057</p>	<p><i>Quick-cfRNA Serum & Plasma Kit</i></p> <p>R1059</p>
DNA/RNA	<p><i>Quick-DNA/RNA Kits</i></p> <p>D7001 (Column) R2130 (Magbead)</p>	<p><i>Quick-cfDNA/cfRNA Serum & Plasma Kit</i></p> <p>R1072</p>

For technical inquiries please contact tech@zymoresearch.com