


Oligo Clean & Concentrator™ Kits

- **Fast & Simple:** Clean and concentrate DNA and RNA oligonucleotides in 2 minutes.
- **Ultra-Pure:** Ready for Next-Gen Sequencing, PCR, ligations, etc.
- **Complete removal:** Eliminate dyes, salts, enzymes, and short oligos.

Improved Recovery of Oligonucleotides



The Oligo Clean & Concentrator™ facilitates >90% recovery of ssDNA oligonucleotides **(A)** and efficient short oligo and nucleotide removal **(B)**.

Product	Cat. No.	Size
Oligo Clean & Concentrator™	D4060	50 preps
	D4061	200 preps
ZR-96 Oligo Clean & Concentrator™	D4062	2 x 96 preps
	D4063	4 x 96 preps