

## Electronic Volume-weight Apparatus



**Model: BL-EVWA-1000**

### **Feature :-**

Grain Volume-Weight Meter is a specialized instrument to measure the grade of food. It can not only measure millet, wheat, sorghum and other small particles of grain, but also corn, soybeans and other large particles of grain.

### **Technical specification :-**

1. Maximum weighing load : 1000g  $\pm$ 2g
2. Weighing resolution : 1g
3. Density barrel volume : 1000 $\pm$ 1.5ml
4. Power : AC 220V $\pm$ 10V 50 Hz
5. Working temperature : 5°C-40°C
7. Relative humidity : <90%RH
6. Blanking mouth diameter : particles 30mm  
large particles 50mm
7. Measuring mode : combined measurement
8. Overall dimension : 530\*230\*405mm
9. Weight : 12g

**SPECIAL  
OFFER**

**75,000 /-**