

● BL-IMS-85 Flake Ice Machine

Description

Biolinkk ice flaking machine BL-IMS-85 is used to produce ice flakes of capacity 85kg/24hr. Features a slide away access door for easy user access to the produced ice. High quality stainless steel body ensures durability.

Features

- * Flake Ice
- * Air Cooling
- * Stainless Steel 304 Exterior
- * Adjustable Feet
- * Slide Away Access Door
- * World Famous Brand Compressor
- * Fully Automatic Controller
- * Improved Design For Higher Reliability



Specifications

Model	BL-IMS-85
Ice Making Capacity	85kg/24h
Power	480W/220V
Storage Capacity	25kg
Product Dimension	400×510×825mm
Water Supplier	Tap Water
Cooling Type	Fan Cooling
Refrigerant	R134a
Ice Shape	Irregular Crushed Ice
Certification	CE