

# ● BL-IMS-20 Flake Ice Machine

## Description

Biolinkk ice flaking machine BL-IMS-20 is used to produce ice flakes of capacity 20kg/24hr. Features a slide away access door for easy user access to the produced ice. High quality stainless steel body ensures durability.

## Features

- \* Flake Ice
- \* Air Cooling
- \* Stainless Steel 304 Exterior
- \* Adjustable Feet
- \* Slide Away Access Door
- \* World Famous Brand Compressor
- \* Fully Automatic Controller
- \* Improved Design For Higher Reliability



## Specifications

Model	BL-IMS-20
Ice Making Capacity	20kg/24h
Power	280W/220V
Storage Capacity	10kg
Product Dimension	330×470×605mm
Water Supplier	Tap Water
Cooling Type	Fan Cooling
Refrigerant	R134a
Ice Shape	Irregular Crushed Ice
Certification	CE