



Literature-References citing the use of Milenia HybriDetect:

Date: July 16th, 2020

RPA

1. Piepenburg O., Williams C.H., Stemple D.L., Armes N.A.
DNA detection using recombination proteins. *PLoS Biol* (2006; 4(7): e204.)
DOI: <https://doi.org/10.1371/journal.pbio.0040204>
2. Rohrman B., Richards-Kortum R., A paper and plastic device for performing recombinase polymerase amplification of HIV DNA. *Lab on a Chip* (2012)
DOI: <https://doi.org/10.1039/c2lc40423k>
3. Crannell Z., Castellanos-Gonzalez A., Irani A., et. al., Nucleic acid test to diagnose cryptosporidiosis: Lab assessment in animal and patient specimens. *Analytical Chemistry* (2014)
DOI: <https://doi.org/10.1021/ac403750z>
4. Cordray MS et al., A Paper and Plastic Device for the Combined Isothermal Amplification and Lateral Flow Detection of Plasmodium DNA. *Malar. J.* (2015; 14, 472)
DOI: <https://doi.org/10.1186/s12936-015-0995-6>
5. Rosser et al., Isothermal Recombinase Polymerase Amplification (RPA) of Schistosoma haematobium DNA and Oligochromatographic Lateral Flow Detection. *Parasit Vectors* (2015; 8)
DOI: <https://doi.org/10.1186/s13071-015-1055-3>
6. Castellanos-Gonzalez A et al., A Novel Test to Diagnose Canine Visceral Leishmaniasis at the Point of Care. *Am. J. Med Hyg.* (2015, 93 (5), 970-975)
DOI: <https://doi.org/10.4269/ajtmh.15-0145>
7. Chao CC et al., Development of Recombinase Polymerase Amplification Assays for Detection of Orienta Tsutsugamushi or Rickettsia typhi. *PLoS Negl. Trop. Dis.* (2015; 9 (7))
DOI: <https://doi.org/10.1371/journal.pntd.0003884>
8. Nair G. et al., Detection of Entamoeba histolytica by Recombinase Polymerase Amplification. *Am. J. Trop. Med.* (2015; 93 (3) 591-595)
DOI: <https://doi.org/10.4269/ajtmh.15-0276>
9. Szanto-Egesz R., Janosi A., Mohr A., et. al., Breed-Specific Detection of Mangalica Meat in Food Products. *Food Analytical Methods* (2016)
DOI: <https://doi.org/10.1007/s12161-015-0261-0>



10. Sun K., Ying W., Yu X., et. al., Recombinase polymerase amplification combined with a lateral flow dipstick for rapid and visual detection of *Schistosoma japonicum*. *Parasites and Vectors* (2016)
DOI: <https://doi.org/10.1186/s13071-016-1745-5>
11. Rohrman B. et al., Inhibition of Recombinase Polymerase Amplification by Background DNA: A Lateral Flow-Based Method for Enriching Target DNA. *Anal. Chem.* (2015 87 (3), 1963-1967)
DOI: <http://pubs.acs.org/doi/pdf/10.1021/ac504365v>
12. Saldarriaga O. A. et al., An Innovative Field-Applicable Molecular Test to Diagnose Cutaneous *Leishmania Viannii* ssp. *Infections. PLoS Negl. Trop. Dis.* (2016; 10 (4))
DOI: <https://doi.org/10.1371/journal.pntd.0004638>
13. Yang Y et al., Development of an Isothermal Amplification-Based Assay for Rapid Visual Detection of an Orf Virus. *Virology J.* (2016; 13:46)
DOI: <https://doi.org/10.1186/s12985-016-0502-x>
14. Hou P., Wang H., Zhao G., et. al., Rapid detection of infectious bovine Rhinotracheitis virus using recombinase polymerase amplification assays. *BMC Veterinary Research* (2017)
DOI: <https://doi.org/10.1186/s12917-017-1284-0>
15. Yang Y., Qin X., Zhang X., et. al., Development of real-time and lateral flow dipstick recombinase polymerase amplification assays for rapid detection of goatpox virus and sheeppox virus. *Virology Journal* (2017)
DOI: <https://doi.org/10.1186/s12985-017-0792-7>
16. Yang Y., Qin X., Sun Y., et. al., Development of Isothermal Recombinase Polymerase Amplification Assay for Rapid Detection of Porcine Circovirus Type 2. *BioMed Research International* (2017)
DOI: <https://doi.org/10.1155/2017/8403642>
17. Liu Y., Wang X., Wei X., et. al., Rapid Authentication of Ginkgo biloba Herbal Products Using the Recombinase Polymerase Amplification Assay. *Scientific Reports* (2018)
DOI: <https://doi.org/10.1038/s41598-018-26402-8>
18. Tu P., A Recombinase Polymerase Amplification Lateral Flow Dipstick for Field Diagnosis of Bovine Leukemia Virus Infection and its Effectiveness Compared to iPCR. and ELISA. *Journal of Antivirals & Antiretrovirus* (2018)
DOI: <https://doi.org/10.4172/1948-5964.1000178>
19. Hou P., Zhao G., Wang H., et. al., Rapid detection of bovine viral diarrhea virus using recombinase polymerase amplification combined with lateral flow dipstick assays in bulk milk. *Veterinarski arhiv* (2018)
DOI: <https://doi.org/10.24099/vet.arhiv.0145>
20. Qi Y., Shao Y., Rao J., et. al., Development of a rapid and visual detection method for *Rickettsia rickettsii* combining recombinase polymerase assay with lateral flow test. *PLoS ONE* (2018)
DOI: <https://doi.org/10.1371/journal.pone.0207811>
21. Zhao G et al. Development of a recombinase polymerase amplification combined with a lateral flow dipstick assay for rapid detection of the *Mycoplasma bovis*. *BMC Vet Res.* (2018) Dec 20;14(1):412
DOI: <https://doi.org/10.1186/s12917-018-1703-x>
22. Zhao G. et al., Rapid visual detection of *Mycobacterium avium* subsp. *paratuberculosis* by recombinase polymerase amplification combined with a lateral flow dipstick. *J Vet Sci.* (2018) Mar 31;19(2):242-250.
DOI: <https://doi.org/10.4142/jvs.2018.19.2.242>



23. Qi Y. et al., Development of a rapid and visual nucleotide detection method for a Chinese epidemic strain of *Orientia tsutsugamushi* based on recombinase polymerase amplification assay and lateral flow test. *Int J Infect Dis.* 2018 May;70:42-50.
DOI: <https://doi.org/10.1016/j.ijid.2018.03.003>
24. Ahmed F. A. et al., Genome-informed diagnostics for specific and rapid detection of *Pectobacterium* species using recombinase polymerase amplification coupled with a lateral flow device. *Sci Rep.* (2018) Oct 29;8(1):15972
DOI: <https://doi.org/10.1038/s41598-018-34275-0>
25. Qi Y. et al., Development of a rapid and visual detection method for *Rickettsia rickettsii* combining recombinase polymerase assay with lateral flow test. *PLoS One.* (2018) Nov 26;13(11)
DOI: <https://doi.org/10.1371/journal.pone.0207811>
26. Wang H. et al., Rapid visual detection of cyprinid herpesvirus 2 by recombinase polymerase amplification combined with a lateral flow dipstick. *J Fish Dis.* (2018) May 28
DOI: <https://doi.org/10.1111/jfd.12808>
27. Wang H. et al., Development and evaluation of serotype-specific recombinase polymerase amplification combined with lateral flow dipstick assays for the diagnosis of foot-and-mouth disease virus serotype A, O and Asia1. *BMC Vet Res.* (2018) Nov 20;14(1)
DOI: <https://doi.org/10.1186/s12917-018-1644-4>
28. Soliman H. et al., Recombinase polymerase amplification assay combined with a lateral flow dipstick for rapid detection of *Tetracapsuloides bryosalmonae*, the causative agent of proliferative kidney disease in salmonids. *Parasit Vectors.* (2018) Apr 11;11(1):234.
DOI: <https://doi.org/10.1186/s13071-018-2825-5>
29. Gootenberg J., Abudayyeh O., Kellner M., et. al., Multiplexed and portable nucleic acid detection platform with Cas13, Cas12a, and Csm6. *Science* (2018)
DOI: [science.sciencemag.org/](https://doi.org/10.1126/science.1262016)
CRISPR
30. James A., Todd S., Pollak N., et. al., Ebolavirus diagnosis made simple, comparable and faster than molecular detection methods: Preparing for the future. *Virology Journals* (2018)
DOI: <https://doi.org/10.1186/s12985-018-0985-8>
31. Myhrvold C., Freije C., Gootenberg J., et. al., Field-deployable viral diagnostics using CRISPR-Cas13. *Science* (2018)
DOI: [science.sciencemag.org/](https://doi.org/10.1126/science.1262016)
CRISPR
32. Zingg J., Daunert S., Trinucleotide Rolling Circle Amplification: A Novel Method for the Detection of RNA and DNA. *Methods and Protocols* (2018)
DOI: <https://doi.org/10.3390/mps1020015>
33. Soliman H. and El-Matoubli M., Rapid detection and differentiation of carp oedema virus and cyprinid herpes virus-3 in koi and common carp. *J Fish Dis.* (2018) May;41(5):761-772.
DOI: <https://doi.org/10.1111/jfd.12774>
34. Kuşkuçcu M.A., Aygün G., Karakullukçu A., Imamova N., Küçükbasmaçlı Ö., Midilli K., Development of a rapid molecular test format based on isothermal recombinase polymerase amplification technique and lateral flow method for detection of bla_{oxa-48}. *Klimik Dergisi* (2018); V 31 (3, 185-189)
DOI: <https://doi.org/10.5152/kd.2018.46>



35. Zasada A. A. et al., Isothermal DNA amplification combined with lateral flow dipsticks for detection of biothreat agents. *Anal Biochem.* (2018) Nov 1;560:60-66.
DOI: <https://doi.org/10.1016/j.ab.2018.09.008>
36. Zhao G. et al., Use of a recombinase polymerase amplification commercial kit for rapid visual detection of *Pasteurella multocida*. *BMC Vet. Res.* 2019; 15:154
DOI: <https://doi.org/10.1371/journal.pbio.0040204>
37. Li T.-T. et al., Rapid and Visual Detection of *Trichinella* Spp. Using a Lateral Flow Strip-Based Recombinase Polymerase Amplification (LF-RPA) Assay. *Front. Cell. Infect. Microbiol.* (2019; 9 1)
DOI: <https://doi.org/10.3389/fcimb.2019.00001>
38. Ma Q. et al., Development of a lateral flow recombinase polymerase amplification assay for rapid and visual detection of *Cryptococcus neoformans*/C. *gattii* in cerebral spinal fluid. *BMC Infect. Dis.* (2019; 19:108)
DOI: <https://doi.org/10.1186/s12879-019-3744-6>
39. Peng Y. et al., Rapid detection of *Burkholderia pseudomallei* with a lateral flow recombinase polymerase amplification assay. *PLoS ONE* 2019; 14(7)
DOI: <https://doi.org/10.1371/journal.pone.0213416>
40. Sullivan, T. J., Dhar, A. K., Cruz-Flores, R., & Bodnar, A. G. Rapid, CRISPR-Based, Field-Deployable Detection Of White Spot Syndrome Virus In Shrimp. *Scientific Reports*, (2019) 9(1),
DOI: <https://doi.org/10.1038/s41598-019-56170-y>
CRISPR
41. Bai J., Lin H., Li H., Zhou Y., Liu J., Zhong G., Cas12a-Based On-Site and Rapid Nucleic Acid Detection of African Swine Fever. *Frontiers in Microbiology* (2019); V 10(1-9)
DOI: <https://doi.org/10.3389/fmicb.2019.02830>
42. Abudayyeh O. O., Gootenberg J. S., Kellner M. J., Zhang F., Nucleic Acid Detection of Plant Genes Using CRISPR-Cas13. *The CRISPR Journal* (2019), V 2 (3, 165-171)
DOI: <https://doi.org/10.1089/crispr.2019.0011>
CRISPR
43. Kellner M. J., Koob J. G., Gootenberg J. S., Abudayyeh O. O., Zhang F., SHERLOCK: nucleic acid detection with CRISPR nucleases. *Nature Protocols* (2019); V 14 (10, 2986-3012)
DOI: <https://doi.org/10.1038/s41596-019-0210-2>
CRISPR
44. Chia J. Y., Tay Y. P., Poh D., Tay B.H., Koh E., Teo J., Yap H. M., Lim M., Soh L. T., Lee Y.K., Lee C. K. 1 Development and evaluation of a rapid on-site water pathogen detection system for water quality monitoring. *Journal of Chemical Information and Modeling* (2019).
DOI: <https://doi.org/10.1101/731778>
45. Li J., Macdonald J., Von Stetten F., Review: a comprehensive summary of a decade development of the recombinase polymerase amplification. *Analyst* (2019); V 144 (1, 31-67)
DOI: <https://doi.org/10.1039/c8an01621f>
46. Peng Y., Zheng X., Kann B., et. al., Rapid detection of *Burkholderia pseudomallei* with a lateral flow recombinase polymerase amplification assay. *PLoS One* (2019)
DOI: <https://doi.org/10.1371/journal.pone.0213416>



47. Tsou J., Leng Q., Jiang F., A CRISPR Test for Detection of Circulating Nuclei Acids. *Translational Oncology* (2019).
DOI: <https://doi.org/10.1016/j.tranon.2019.08.011>
CRISPR
48. Zhang T., Liu M., Yin R., et. al., Rapid and simple detection of *Glaesserella parasuis* in synovial fluid by recombinase polymerase amplification and lateral flow strip. *BMC Veterinary Research* (2019)
DOI: <https://doi.org/10.1186/s12917-019-2039-x>
49. Dai T., Hu T., Yang X., et. al., A recombinase polymerase amplification lateral flow dipstick assay for rapid detection of the quarantine citrus pathogen in China, *Phytophthora hibernalis*. *PeerJ* (2019)
DOI: <https://doi.org/10.7717/peerj.8083>
50. Dai T., Yang X., Hu T., et. al., Comparative evaluation of a novel recombinase polymerase amplification-lateral flow dipstick (Rpa-lfd) assay, lamp, conventional pcr, and leaf-disc baiting methods for detection of *phytophthora sojae*. *Frontiers in Microbiology* (2019)
DOI: <https://doi.org/10.3389/fmicb.2019.01884>
51. Wu L., Ye L., Wang Z., et. al., Utilization of recombinase polymerase amplification combined with a lateral flow strip for detection of *Perkinsus behaiensis* in the oyster *Crassostrea hongkongensis*. *Parasites & vectors* (2019)
DOI: <https://doi.org/10.1186/s13071-019-3624-3>
52. Tsou J., Leng Q., Jiang F., A CRISPR Test for Detection of Circulating Nuclei Acids. *Translational Oncology* (2019)
DOI: <https://doi.org/10.1016/j.tranon.2019.08.011>
CRISPR
53. Curti L., Pereyra-Bonnet F., Gimenez C., An ultrasensitive, rapid, and portable coronavirus SARS-CoV-2 sequence detection method based on CRISPR-Cas12. *bioRxiv* (2020); V 11 (2020.02.29.971127)
DOI: <https://doi.org/10.1101/2020.02.29.971127>
CRISPR
54. Metsky H. C., Freije C. A., Kosoko-Thoroddsen T-S. F., Sabeti P.C., Myhrvold C., CRISPR-based surveillance for COVID-19 using genomically-comprehensive machine learning design. *bioRxiv* (2020), 2020.02.26.967026
DOI: <https://doi.org/10.1101/2020.02.26.967026>
CRISPR
55. Zhang F., Abudayyeh O. O., Gootenberg J.S., Mathers L. A protocol for detection of COVID-19 using CRISPR diagnostics. *Broad Communications* (2020), [https://www.broadinstitute.org/files/publications/special/COVID-19%20detection%20\(updated\).pdf](https://www.broadinstitute.org/files/publications/special/COVID-19%20detection%20(updated).pdf)
CRISPR
56. Qian J., Boswell S., Chidley C., et. al., An enhanced isothermal amplification assay for viral detection. *BioRxiv* preprint (2020)
DOI: <https://doi.org/10.1101/2020.05.28.118059>
57. Gao X., Liu X., Zhang Y., et. al., Rapid and visual detection of porcine deltacoronavirus by recombinase polymerase amplification combined with a lateral flow dipstick. *BMC Veterinary Research* (2020)
DOI: <https://doi.org/10.1186/s12917-020-02341-3>



58. Li Z., Torres J., Goossens J., et. al., Development of a recombinase polymerase amplification lateral flow assay for the detection of active *Trypanosoma evansi* infections. *PLoS Neglected Tropical Diseases* (2020)
DOI: <https://doi.org/10.1371/journal.pntd.0008044>
59. Dai T., Wang A., Yang X., Yu X., Tain W., Xu Y., PHYCI _ 587572 : An RxLR Effector Gene and New Biomarker in A Recombinase Polymerase Amplification Assay for Rapid Detection of *Phytophthora cinnamomi*. *Forests* (2020); V 11 (306)
DOI: <https://doi.org/10.3390/f11030306>
60. Ravindran V. B., Khallaf B., Surapaneni A., Crosbie N.D., Soni S.K., Ball A.S., Detection of Helminth Ova in Wastewater Using Recombinase Polymerase Amplification Coupled to Lateral Flow Strips. *Water* (2020), V 12 (3, 691)
DOI: <https://doi.org/10.3390/w12030691>
61. Wang X., Xion E., Tian T., Cheng M., Lin W., Wang H. Supporting Information Clustered Regularly Interspaced Short Palindromic Repeats / Cas9- Mediated Lateral Flow Nucleic Acid Assay. *ACS Nano* (2020), IF 13.903
DOI: <https://doi.org/10.1021/acsnano.0c00022>
CRISPR
62. Jiang W., Ren Y., Han X., Xue J., Shan T., Chen Z., Liu Y., Quan W. Recombinase polymerase amplification-lateral flow (RPA-LF) assay combined with immunomagnetic separation for rapid visual detection of *Vibrio parahaemolyticus* in raw oysters. *Analytical and Bioanalytical Chemistry* (2020), V 51
DOI: <https://doi.org/10.1007/s00216-020-02532-9>
63. Azhar M. et al., Rapid, field-deployable nucleobase detection and identification using FnCas9. *bioRxiv* preprint (2020)
DOI: <https://doi.org/10.1101/2020.04.07.028167>
CRISPR
64. Baerwald M., Goodbla A., nagarajan R., Gootenberg J., Abudayyeh O., Zhang F., Schreier A., Rapid and accurate species identification for ecological studies and monitoring using CRISPR-based SHERLOCK. *Molecular Ecology Resources* (2020)
DOI: <https://doi.org/10.1111/1755-0998.13186>
CRISPR
65. Curti L., Pereyra-Bonnet F., Repizo G., et. al., CRISPR-based platform for carbapenemases and emerging viruses detection using Cas12a (Cpf1) effector nuclease. *Emerging Microbes and Infections* (2020), vol 9 (1)
DOI: <https://doi.org/10.1080/22221751.2020.176385>
CRISPR
66. Hao L., Zhao R., Ngambenjawong C., et. al., CRISPR-Cas-amplified urine biomarkers for multiplexed and portable cancer diagnostics, *bioRxiv* preprint (2020)
DOI: <https://doi.org/10.1101/2020.06.17.157180>
CRISPR
67. Hu J., Wang Y., Su H., et. al., Rapid analysis of *Escherichia coli* O157:H7 using isothermal recombinase polymerase amplification combined with triple-labeled nucleotide probes. *Molecular and Cellular Probes* (2020) vol 50.
DOI: <https://doi.org/10.1016/j.mcp.2019.101501>



68. James A., Alwneh J., COVID-19 infection diagnosis: Potential impact of isothermal amplification technology to reduce community transmission of SARS-CoV-2. *Diagnostics* (2020)
DOI: <https://doi.org/10.3390/diagnostics10060399>
69. Joung J., Ladha A., Saito M., et. al., Point-of-care testing for COVID-19 using SHERLOCK diagnostics. *Medrxiv* preprint (2020)
DOI: <https://doi.org/10.1101/2020.05.04.20091231>
70. Ju Y., Li C., Shen P., et. al., Rapid and visual detection of *Verticillium dahliae* using recombinase polymerase amplification combined with lateral flow dipstick. *Crop Protection* (2020)
DOI: <https://doi.org/10.1016/j.cropro.2020.105>
71. Kaminski M., Alcantar M., Lape I., et.al., A CRISPR-based assay for the detection of opportunistic infections post-transplantation and for the monitoring of transplant rejection. *Nature Biomedical Engineering* (2020)
DOI: <https://doi.org/10.1038/s41551-020-0546-5>
CRISPR
72. Kanitchinda S., Srisala J., Suebsing R., et. al., CRISPR-Cas fluorescent cleavage assay coupled with recombinase polymerase amplification for sensitive and specific detection of *Enterocytozoon hepatopenaei*. *Biotechnology Reports* (2020)
DOI: <https://doi.org/10.1016/j.btre.2020.e00485>
CRISPR
73. Lalremruata A., Nguyen T., McCall M., et. al., Recombinase Polymerase Amplification and Lateral Flow Assay for Ultrasensitive Detection of Low-Density. *Journal of Clinical Microbiology* (2020)
DOI: <https://doi.org/10.1128/JCM.01879-19>
74. Lu L., Li F., Chen Q., et. al., Rapid detection of African swine fever virus using Cas12a-based portable paper diagnostics. *Cell Discovery* (2020)
DOI: <https://doi.org/10.1038/s41421-020-0151-5>
CRISPR
75. Aman R., Mahas A., Mahfouz M., Nucleic Acid Detection Using CRISPR/Cas Biosensing Technologies. *ACS synthetic biology* (2020)
DOI: <https://doi.org/10.1021/acssynbio.9b00507>
CRISPR
76. Chang Y., Deng Y., Li T., Visual detection of porcine reproductive and respiratory syndrome virus using CRISPR-Cas13a. *Transboundary and Emerging Diseases* (2020) Vol 67,
DOI: <https://doi.org/10.1111/tbed.13368>
CRISPR
77. Yoshimi K., Takeshita K., Yamayoshi S., et. al., Rapid and accurate detection of novel coronavirus SARS-CoV-2 using CRISPR-Cas3. *medRxiv* preprint (2020)
DOI: <https://doi.org/doi.org/10.1101/2020.06.02.20119875>
CRISPR
78. Xinyu L., Zheng Y., Zhang F., A rapid, equipment-free method for detecting *Phytophthora infestans* in the field using a lateral flow strip-based recombinase polymerase amplification assay. *Plant Disease* (2020)
DOI: <https://doi.org/10.1094/PDIS-01-20-0203-SC>



LAMP

1. Kiatpathomchai W., Jaroenram W., Arunrut N., Jitrapakdee S., Flegel TW. Shrimp Taura syndrome virus detection by reverse transcription loop-mediated isothermal amplification combined with a lateral flow dipstick. *Journal of Virological Methods* (2008; 153: 214–217)
DOI: <https://doi.org/10.1016/j.jviromet.2008.06.025>
2. Kikuchi, T., Aikawa, T., Oeda, Y., Karim, N., & Kanzaki, N. A rapid and precise diagnostic method for detecting the pinewood nematode *Bursaphelenchus xylophilus* by loop-mediated isothermal amplification. *Phytopathology*, 2009, 99(12), 1365–1369.
DOI: <https://doi.org/10.1094/PHYTO-99-12-1365>
3. Soliman H., El-Matbouli M., Loop mediated isothermal amplification combined with nucleic acid lateral flow strip for diagnosis of cyprinid herpes virus-3. *Molecular and Cellular Probes* (2010)
DOI: <https://doi.org/10.1016/j.mcp.2009.09.002>
4. Ding W. C., Chen J., Shi Y. H., Lu X. J., Li M. Y. Rapid and sensitive detection of infectious spleen and kidney necrosis virus by loop-mediated isothermal amplification combined with a lateral flow dipstick. *Arch. Virol.* (2010; 155: 385-389)
DOI: <https://doi.org/10.1007/s00705-010-0593-4>
5. Surasilp T., Longyant S., Rukpratanporn S., Sridulyakul P., Sithigorngul P., Chaivisuthangkura P. Rapid and sensitive detection of *Vibrio vulnificus* by loop-mediated isothermal amplification combined with lateral flow dipstick targeted to rpoS gene. *Molecular and Cellular Probes* (2011; 25: 158-163)
DOI: <https://doi.org/10.1016/j.mcp.2011.04.001>
6. Kaewphinit T., Arunrut N., Kiatpathomchai W., Santiwatanakul S., Jaratsing P., Chansiri K. Detection of *Mycobacterium tuberculosis* by using loop-mediated isothermal amplification combined with a lateral flow dipstick in clinical samples. *Biomed. Res. Int.* (2013; 1-6)
DOI: <https://doi.org/10.1155/2013/926230>
7. Vaagt F et al. Loop-mediated isothermal amplification (LAMP)-based method for rapid mushroom species identification. *J. Agric. Food Chem.* 2013; 61 (8), 1833-1840
DOI: <https://doi.org/10.1021/jf304824b>
8. Sun YL et al., Visual Detection of Canine Parvovirus Based on Loop-Mediated Isothermal Amplification Combined with Enzyme-Linked Immunosorbent Assay and with Lateral Flow Dipstick. *J. Vet. Med. Sci.* (2014)
DOI: <https://doi.org/10.1292/jvms.13-044>
9. Rigano AL et al., Rapid and sensitive detection of *Candidatus Liberibacter asiaticus* by loop mediated isothermal amplification combined with a lateral flow dipstick *BMC Microbiol.* 2014; 1186,
DOI: <http://www.biomedcentral.com/1471-2180/14/86>
10. Kusumawati A. et al., Use of reverse transcription loop-mediated isothermal amplification combined with lateral flow dipstick for an easy and rapid detection of Jembrana disease virus. *Virusdisease.* (2015 Sep;26(3):189-95)
DOI: <https://doi.org/10.1007/s13337-015-0277-5>



11. Yun X., Hao Y., Peng Z., Repaid and accurate detection of *Khuskia oryzae* by loopmediated isothermal amplification combined with a lateral flow dipstick. *Chemical Engineering Transactions* (2015)
DOI: <https://doi.org/10.3303/CET1546147>
12. Yang Y et al., Rapid and sensitive detection of *Babesia bovis* and *Babesia bigemina* by loop-mediated isothermal amplification combined with a lateral flow dipstick, *Vet. Parasitol* (2016; 219, 71-76)
DOI: <https://doi.org/10.1016/j.vetpar.2016.02.004>
13. Huang HL et al., The Development of Loop-Mediated Isothermal Amplification Combined with Lateral Flow Dipstick for Detection of *Karlodinium veneficum*. *Harmful Algae* (2017; 62, 20-29)
DOI: <https://doi.org/10.1016/j.hal.2016.11.022>
14. Kumvongpin R., Jearanaikoon P., Wilailuckana C., et. al., Detection assay for HPV16 and HPV18 by loop-mediated isothermal amplification with lateral flow dipstick tests. *Molecular Medicine Reports* (2017)
DOI: <https://doi.org/10.3892/mmr.2017.6370>
15. Huang H. L. et al., Detection of *Skeletonema costatum* based on loop-mediated isothermal amplification combined with lateral flow dipstick. *Mol Cell Probes*. (2017 Dec;36:36-42)
DOI: <https://doi.org/10.1016/j.mcp.2017.08.003>
16. Lalle M. et al., Loop-Mediated Isothermal Amplification-Lateral-Flow Dipstick (LAMP-LFD) to detect *Toxoplasma gondii* oocyst in ready-to-eat salad. *Food Microbiol*. (2018 Apr;70:137-142)
DOI: <https://doi.org/10.1016/j.fm.2017.10.001>
17. Kanchanaphum P., Time Course of Detection of Human Male DNA from Stained Blood Sample on Various Surfaces by Loop Mediated Isothermal Amplification and Polymerase Chain Reaction. *Biomed Res Int*. (2018 Mar 22;2018:2981862)
DOI: <https://doi.org/10.1155/2018/2981862>
18. Mamba T. S. et al., Lateral Flow Loop-Mediated Isothermal Amplification Test with Stem Primers: Detection of *Cryptosporidium* Species in Kenyan Children Presenting with Diarrhea. *J Trop Med*. (2018 Feb 26;2018:7659730)
DOI: <https://doi.org/10.1155/2018/7659730>
19. Zasada A. A. et al., Isothermal DNA amplification combined with lateral flow dipsticks for detection of biothreat agents. *Anal Biochem*. (2018 Nov 1;560:60-66)
DOI: <https://doi.org/10.1016/j.ab.2018.09.008>
20. Cai C., Jung C., Bhadra S., Ellington A.D., Phosphorothioated Primers Lead to Loop-Mediated Isothermal Amplification at Low Temperatures. *Analytical Chemistry* (2018, 90, 82900-8294)
DOI: <https://pubs.acs.org/doi/10.1021/acs.analchem.8b02062>
21. Zhai X., A loop-mediated isothermal amplification coupling with a lateral flow dipstick for rapid and specific detection of fowl adenovirus serotype-4. *J. Virol. Methods* (2019; 270, 79-86)
DOI: <https://doi.org/10.1016/j.jviromet.2019.04.026>
22. Sagcan H., Turgut Kara N., Detection of Potato ring rot Pathogen *Clavibacter michiganensis* subsp. *sepedonicus* by Loop-mediated isothermal amplification (LAMP) assay. *Scientific Reports* (2019)
DOI: <https://doi.org/10.1038/s41598-019-56680-9>



23. Thongphueak D. et al., Development of the Rapid Test Kit for the Identification of *Campylobacter* spp. Based on Loop-mediated Isothermal Amplification (LAMP) in Combination with a Lateral Flow Dipstick (LFD) and Gold Nano-DNA Probe (AuNPs). *Sci&Tech. Asia* (2019; 24(1) 63-71)
DOI: <https://doi.org/10.14456/scitechasia.2019.7>
24. Ledlod, S., Bunroddith, K., Areekit, S., Santiwatanakul, S., & Chansiri, K. (2019). Development of a duplex lateral flow dipstick test for the detection and differentiation of *Listeria* spp. and *Listeria monocytogenes* in meat products based on loop-mediated isothermal amplification. *Journal of Chromatography B. Analytical Technologies in the Biomedical and Life Sciences*, 1139(October 2019), 121834.
DOI: <https://doi.org/10.1016/j.jchromb.2019.121834>
25. Bai, J., Lin, H., Li, H., Zhou, Y., Liu, J., Zhong, G., ... Huang, L. Cas12a-Based On-Site and Rapid Nucleic Acid Detection of African Swine Fever. *Frontiers in Microbiology*, 2019, 10(December), 1–9.
DOI: <https://doi.org/10.3389/fmicb.2019.02830>
CRISPR
26. Abudayyeh, O. O., Gootenberg, J. S., Kellner, M. J., & Zhang, F. Nucleic Acid Detection of Plant Genes Using CRISPR-Cas13. *The CRISPR Journal*, (2019) 2(3), 165–171.
DOI: <https://doi.org/10.1089/crispr.2019.0011>
CRISPR
27. Kellner, M. J., Koob, J. G., Gootenberg, J. S., Abudayyeh, O. O., & Zhang, F. SHERLOCK: nucleic acid detection with CRISPR nucleases. *Nature Protocols*, (2019) 14(10), 2986–3012.
DOI: <https://doi.org/10.1038/s41596-019-0210-2>
CRISPR
28. Huang H. L., Gao W.F., Zhu P., Zhou C.X., Qiao L.L., Dang C.Y., Pang J.H., Yan X.J., Molecular method for rapid detection of the red tide dinoflagellate *Karenia mikimotoi* in the coastal region of Xiangshan Bay, China. *Journal of Microbiological Methods* (2020); V 168 (105801)
DOI: <https://doi.org/10.1016/j.mimet.2019.105801>
29. Nguyen L., Smith B., Jain P., Enhancement of trans-cleavage activity of Cas12a with engineered crRNA enables amplified nucleic acid detection. *BioRxiv* preprint, (2020).
DOI: <https://doi.org/10.1101/2020.04.13.036079>
CRISPR
30. Broughton J. P., Deng X., Fasching C.L., Singh J., Streithorst J., Granados, A., Sotomayor-Gonzalez A., Zorn K., Gopez A., Hsu E., Gu W., Miller S., Pan C., Guevara H., Wadford D., Chen J., Chio C.Y., Rapid Detection of 2019 Novel Coronavirus SARS-CoV-2 Using a CRISPR-based DETECTR Lateral Flow Assay. *medRxiv* (2020); V 415 (2020.03.06.20032334)
DOI: <https://doi.org/10.1101/2020.03.06.20032334>
CRISPR
31. Broughton J. P., Wayne D., Fasching C. L., Singh J., Charles Y., Chen J. S., A protocol for rapid detection of the 2019 novel coronavirus SARS-CoV-2 using CRISPR diagnostics : SARS-CoV-2 DETECTR. *MammothBiosciences* (2020)
<https://mammoth.bio/2020/02/15/white-paper-a-protocol-for-rapid-detection-of-sars-cov-2-using-crispr-sars-cov-2-detectr/>
CRISPR



32. Ađel H.E., Sađcan H., Ceyhan I., Durmaz R., Optimization of the Isothermal Amplification Method for Mycobacterium Tuberculosis Detection and Visualization for Fieldwork. *Turk J Med Sci* (2020)
DOI: <https://doi.org/10.3906/sag-1910-6>
33. Zhang Y., Zhang Y., Xie K., Evaluation of CRISPR/Cas12a-based DNA detection for fast pathogen diagnosis and GMO test in rice. *Molecular Breeding* (2020); V 40,11
DOI: <https://doi.org/10.1007/s11032-019-1092-2>
CRISPR
34. Li H., Yuan G., Luo Y., Yu Y., Ai T., Su J., Rapid and sensitive detection of infectious spleen and kidney necrosis virus by recombinase polymerase amplification combined with lateral flow dipsticks. *Aquaculture* (2020); V 519, 734926
DOI: <https://doi.org/10.1016/j.aquaculture.2020.734926>
35. Yang X., Xie J., Hu S., Zhan W., Duan L., Chen K., Zhang C., Yin A., Luo M., Rapid and visual detection of enterovirus using recombinase polymerase amplification combined with lateral flow strips. *Sensors and Actuators B: Chemical* (2020); V311, 127903
DOI: <https://doi.org/10.1016/j.snb.2020.127903>
36. Qian J. et al, An enhanced isothermal amplification assay for viral detection. *bioRxiv*, 2020
DOI: <https://doi.org/10.1101/2020.05.28.118059>
37. Allgöwer S., Hartmann C., Holzhauser T., The development of highly specific and sensitive primers for the detection of potentially allergenic soybean (*Glycine max*) using loop-mediated isothermal amplification combined with lateral flow dipstick (LAMP-LFD). *Foods* (2020) vol: 9 (4)
DOI: <https://doi.org/10.3390/foods9040423>
38. Brogton J.P., et al., CRISPR-Cas12-based detection of SARS-CoV-2. *Nature biotechnology* (2020)
DOI: <https://doi.org/10.1038/s41587-020-0513-4>
CRISPR
39. Santiago I., Trends and Innovations in Biosensors for COVID-19 Mass Testing, *ChemBioChem* (2020)
DOI: <https://doi.org/10.1002/cbic.202000250>
40. Nguyen L., Smith B., Jain P., Enhancement of trans-cleavage activity of Cas12a with engineered crRNA enables amplified nucleic acid detection. *bioRxiv* preprint (2020)
DOI: <https://doi.org/10.1101/2020.04.13.036079>
CRISPR
41. Schermer B., Fabretti F., Damagnez M., et. al., Rapid SARS-CoV-2 testing in primary material based on a novel multiplex LAMP assay, *medRxiv* preprint (2020)
DOI: <https://doi.org/10.1101/2020.06.18.20130377>
42. Zasada A., Wiatrzyk A., Czajka U., et. Al., Application of loop-mediated isothermal amplification combined with colorimetric and lateral flow dipstick visualization as the potential point-of-care testing for *Corynebacterium diphtheria*. *BMC Infectious Diseases* (2020 vol 20(1))
DOI: <https://doi.org/10.1186/s12879-020-05037-z>

PCR

1. Wong YL et al., Rapid authentication of Cordyceps by lateral flow dipstick. *J. Pharm. Biomed. Anal* (2015, 111, 306-310)
DOI: <https://doi.org/10.1016/j.jpba.2015.04.003>
2. Muangszwan W., Ruangsuj P., Chaichanachaicharn P., et. al., A novel nucleic lateral flow assay for screening of PHA-producing haloarchaea. *Journal of Microbiological Methods* (2019)
DOI: <https://doi.org/10.1016/j.mimet.2015.06.012>
3. Taboada L. et al., A new method for the rapid detection of Atlantic cod (*Gadus morhua*), Pacific cod (*Gadus macrocephalus*), Alaska pollock (*Gadus chalcogrammus*) and ling (*Molva molva*) using a lateral flow dipstick assay. *Food Chem.* (2017; 233, 182-189)
DOI: <https://doi.org/10.1016/j.foodchem.2017.04.087>
4. Fischer C., Wessels H., Paschke-Kratzin A., et. al., Aptamers: Universal capture units for lateral flow applications. *Analytical Biochemistry* (2017).
DOI: <https://doi.org/10.1016/j.ab.2017.01.012>
5. Breitbach A., Jacob F., Hutzler M., Koob J., Simple Rapid Test for Establishing the Presence of Beer-Spoilage Organisms. *Brauwelt International* (2017)
6. Pecchia S. and Da Lio D., Development of a rapid PCR-Nucleic Acid Lateral Flow Immunoassay (PCR-NALFIA) based on rDNA IGS sequence analysis for the detection of *Macrophomina phaseolina* in soil. *Journal of Microbiological Methods* (2018, Vol 151)
DOI: <https://doi.org/10.1016/j.mimet.2018.06.010>
7. Janz V., Schoon J., Morgenstern C., Rapid detection of periprosthetic joint infection using a combination of 16s rDNA polymerase chain reaction and lateral flow immunoassay. *Bone Joint Research* (2018)
DOI: <https://doi.org/10.1302/2046-3758.71.BJR-2017-0103.R2>
8. Zaky, W. I., Tomaino, F. R., Pilotte, N., Laney, S. J., & Williams, S. A. Backpack PCR: A point-of-collection diagnostic platform for the rapid detection of *Brugia* parasites in mosquitoes. *PLoS Neglected Tropical Diseases*, (2018, 12(11), 1–17)
DOI: <https://doi.org/10.1371/journal.pntd.0006962>
9. Seo H., Kil E., Fadhila C., et. al., Rapid diagnosis of two marine viruses, red sea bream iridovirus and viral hemorrhagic septicemia virus by PCR combined with lateral flow assay. *VirusDisease* (2020)
DOI: <https://doi.org/10.1007/s13337-020-00577-z>