



PRODUCT INFORMATION

PSB Powder (PD0100)

Product name: PBS Powder (Phosphate buffered saline),

Catalog Number: PD0100

Product Description

This PBS formulation is recommended for use in immunohistology protocols.

Each package makes 1 L of 10 concentrated PBS or 10 liters of 1 x solution. Each 1 x PBS solution contains 137 mM Sodium Chloride, 2.7mM Potassium Chloride and 10mM Phosphate Buffer, pH 7.3-7.5(1.0%, Water) at 25°C.

Storage Temperature: Room temperature.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.