



PRODUCT INFORMATION

PBS Blocking Buffer with BSA Premixed Powder

Product information for SB0626:

Description:

BBI PBS Blocking Buffer with BSA Premixed Powder is used for immunodetection assays such as Western blots, dot blots and ELISA.

1 pack will make 1 liter of 1X strength blocking solution contain 3% BSA.

Compositions:

SB0626 contains:

Na₂HPO₄·7H₂O

KH₂PO₄

NaCl

KCl

BSA

Instructions:

1. Dissolve 1 pack of powder (~40.6g) in 1 liter double-distilled water;
2. Used as blocking buffer in experiments, e.g. western blot, ELISA.

Storage condition:

The reagent is shipped at ambient temperature. Upon arrival, store the reagent at 4°C. When stored and used properly, the reagent is stable for 1 year.

Note:

For research, laboratory or further evaluation purposes. NOT FOR HUMAN USE.
