

High-voltage Power Supply (BL-1000E)



| Technical specifications | |
|--------------------------|------------------------------------|
| Cat. No. | BL-1000E |
| Type of Output | Constant-Voltage, Current or Power |
| Increment | 1V, 1mA, 1W |
| Output Range | 30--1000V, 1--500mA, 1--300W. |
| Timer Range | 1min~99hr, 59min |
| V-hr Range | 1~99999V-hr/ Increment 1V-hr |
| Display | LCD display with backlight |
| Output jacks | Four sets of output jacks |

Characteristic features:

- Touch keys, dual core microprocessor intelligent control;
- 0-9 digital button, make the operation more simple and convenient;
- Indicate the preset value and the actual output value at the same time;
- It can store 100 electrophoresis methods;
- With automatic memory function;
- With step-by-step control function (editable 10 groups, each group of up to make 10 program automatic connection operation);
- With a pause control function (including real-time tuning)
- With standard, timing, V-hr, step-by-step operation functions;
- With constant voltage, constant current, constant power, disoperation, fault intelligent prompt functions;
- Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating.