

Basic Power Supply (BL-300)



| Technical specifications | |
|--------------------------|------------------------------------|
| Cat. No. | BL-300 |
| Type of Output | Constant-Voltage, Current or Power |
| Increment | 1V, 1mA |
| Output Range | 5~300V, 1~300mA, maximum 90W |
| Timer Range | 1min~99hr, 59min |
| V-hr Range | 1~99999V-hr/ Increment 1V-hr |
| Display | LED |
| Output jacks | Four sets of output jacks |

Characteristic features:

- Molding shell, touch keys, microprocessor intelligent control;
- With automatic memory function;
- With standard, timing operation function;
- Automatic detection of no-load, over-load.