

## Basic Power Supply (BL-600)



| Technical specifications |                                    |
|--------------------------|------------------------------------|
| Cat. No.                 | BL-600                             |
| Type of Output           | Constant-Voltage, Current or Power |
| Increment                | 1V, 1mA                            |
| Output Range             | 5~600V, 1~200mA, maximum 120W      |
| Timer Range              | 1min~99hr, 59min                   |
| V-hr Range               | 1~99999V-hr/ Increment 1V-hr       |
| Display                  | LED                                |
| Output jacks             | Four sets of output jacks          |

### Characteristic features:

- Molding shell, touch keys, microprocessor intelligent control;
- With automatic memory function;
- With standard, timing operation function;
- Automatic detection of no-load, over-load.