

Universal Power Supply (BL-600HC)



Technical specifications	
Cat. No.	BL-600HC
Type of Output	Constant-Voltage, Current or Power
Increment	1V, 1mA, 1W
Output Range	10--600V, 1--500mA, 1--300W.
Timer Range	1min~99hr, 59min
V-hr Range	1~99999V-hr/ Increment 1V-hr
Display	LCD display with backlight
Output jacks	Four sets of output jacks

Characteristic features:

- Molding shell, touch keys, dual core microprocessor intelligent control;
- Indicate the preset value and the actual output value at the same time;
- It can store 10 electrophoresis methods;
- With automatic memory function;
- With standard, timing operation function;
- With constant voltage, constant current, constant power, disoperation, fault intelligent prompt functions;
- Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating.