

High-voltage Power Supply (BL- ECPT3000)



Technical specifications	
Cat. No.	BL- ECPT3000
Type of Output	Constant-Voltage, Current or Power
Increment	1V, 1mA, 1W
Output Range	20--3000V, 1--200mA, 1--200W.
Timer Range	1min~99hr, 59min
V-hr Range	1~99999V-hr/ Increment 1V-hr
Display	LCD display with backlight
Output jacks	Four sets of output jacks

Characteristic features:

- Touch keys, dual core microprocessor intelligent control;
- Indicate the preset value and the actual output value at the same time;
- It can store 10 electrophoresis methods;
- With the temperature control operation, range of temperature is RT+10°C~65°C;
- With automatic memory function;
- With Standard, Timing, Temperature control operation function;
- With constant voltage, constant current, constant power, disoperation, fault intelligent prompt functions;
- Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating.