

## Horizontal Electrophoresis (BL- SP7)



### Technical specifications:

Gel Size (W×L)	70×250mm×2; 90×250mm×2.
Buffer Volume	1000ml
Dimension (L×W×H)	350×260×65 (mm)
Net weight	1.6 (kg)

### Characteristic features:

- The BL-SP7 is updated from the BL -SP7A with a double-bridge style design;
- With a cooling device, which can make the electrophoresis process under constant temperature.
- You can adjusting the position of the adjusting rod according to the medium length, to meet the different electrophoretic distance.
- Unfolded electrodes make repair and parts replacement easy and quick.
- Auto-switch-off when the lid is removed.