


Horizontal Electrophoresis (BL- SP7A)



Technical specifications:

Gel Size (W×L)	66×250mm; 82×250mm
Buffer Volume	560ml
Dimension (L×W×H)	318×175×65 (mm)
Net weight	1.0 (kg)

Characteristic features:

- With a cooling device, which can make the electrophoresis process under constant temperature;
- You can adjusting the position of the adjusting rod according to the medium length, to meet the different electrophoretic distance.
- Unfolded electrodes make repair and parts replacement easy and quick.
- Auto-switch-off when the lid is removed.