


Horizontal Electrophoresis (BL- SPAT)



Technical specifications:

Gel Size (W×L)	48×75mm
Sample volume	3, 5, 9.
Buffer Volume	150ml
Dimension (L×W×H)	250×90×50 (mm)
Net weight	0.4 (kg)

Characteristic features:

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make repair and parts replacement easy and quick.
- With the sample background color design on the tray, convenient to adding the sample.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.