


Horizontal Electrophoresis (BL- SPBT)



Technical specifications:

Gel Size (W×L)	78×100mm
Sample volume	10, 15
Buffer Volume	400ml
Dimension (L×W×H)	365×105×55 (mm)
Net weight	0.6 (kg)

Characteristic features:

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make repair and parts replacement easy and quick.
- With the sample background color design on the tray, convenient to adding the sample.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.