


## Horizontal Electrophoresis (BL- SPFT)



### Technical specifications:

Gel Size (W×L)	130×130mm; 78×100mm; 48×75mm.
Sample volume	3, 5, 10, 15, 13, 26.
Buffer Volume	600ml
Dimension (L×W×H)	305×152×80 (mm)
Net weight	1.7 (kg)

### Characteristic features:

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make repair and parts replacement easy and quick.
- With the fluorescent ruler on the tray, can be observed directly into the ultraviolet.
- With a special gel casting, preventing the gel leakage, operation is simple and convenient.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.