


Horizontal Electrophoresis (BL- SPGT)



Technical specifications:

Gel Size (W×L)	200×200mm; 200×150mm; 200×100mm.
Sample volume	17, 22, 36, 44.
Buffer Volume	1800ml
Dimension (L×W×H)	397×230×93 (mm)
Net weight	3.3 (kg)

Characteristic features:

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make repair and parts replacement easy and quick.
- With a special gel casting, preventing the gel leakage, operation is simple and convenient.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.