


Horizontal Electrophoresis (BL- SPHT)



Technical specifications:

Gel Size (W×L)	245×250mm.
Sample volume	19~364
Buffer Volume	2000ml
Dimension (L×W×H)	487×300×93 (mm)
Net weight	3.8 (kg)

Characteristic features :

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Removable electrode enables maintain and replacement easy and quick.
- Adding sample more quick and suitable for the multi-channel pipettes.
- With a special gel casting, preventing the gel leakage, operation is simple and convenient.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.