

## Tank transfer (BL- ZY5)



### Technical specifications

Trans area (W×L)	100×100 mm
Buffer Volume	750 (ml)
Dimension (L×W×H)	150×120×115 (mm)
Net weight	0.8 (kg)

### Characteristic features:

- Electrode spacing is small, greatly save transfer time.
- With the color marking on the splint, to ensure put the gel right into the inner core.
- Unfolded electrodes make repair and parts replacement easy and quick.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.
- Not only can be used as a device independent, and can be used as a module with the buffer tank of BL-SCZ2+.