


Tank transfer (BL- ZY7)



Technical specifications:

Trans area (W×L)	216×110 mm
Buffer Volume	1800 (ml)
Dimension (L×W×H)	290×170×155 (mm)
Net weight	1.7 (kg)

Characteristic features:

- With the color marking on the splint, to ensure put the gel right into the inner core.
- Unfolded electrodes make repair and parts replacement easy and quick.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.
- Not only can be used as a device independent, and can be used as a module with the buffer tank of BL-JX5 and BL-CZ-B.