

High Current Power Supply

CE



Voltage	Current	Power
300 V	2000 mA	200 W
Application Recommended:		
Mini vertical electrophoresis		
Horizontal electrophoresis		
Mini tank transfer (normal/fast)		
Large tank transfer		
Semi-Dry electrophoretic transfer		

Specifications:

- Output voltage: 10 – 300 V in 1v steps
- Output current: 10 – 2000 mA in 10 mA steps
- Output accuracy: $\pm 1\%$
- Output Power: 200W max
- Timing: 1 – 999 min fully adjustable
- Display: 3-digital LED
- Type of output: Constant voltage or Constant current
- Resume function: Yes, after power failure
- Operation conditions: 0-40°C, 0-95% humidity
- Safety features: No load detection, sudden load change detection, overload / short circuit protection.
- Dimension: 31×24×8 cm (L×W×H)
- Weight: 1.5kg
- Certificate: CE

High Voltage Power Supply

CE



Voltage	Current	Power
600 V	800 mA	200 W
Application Recommended:		
Large vertical electrophoresis		
Mini vertical electrophoresis		
Horizontal electrophoresis		
Mini tank transfer (normal)		

Specifications:

- Output voltage: 10 – 600 V in 1v steps
- Output current: 10 – 800 mA in 1 mA steps
- Output accuracy: $\pm 1\%$
- Output Power: 200W max
- Timing: 1 – 999 min fully adjustable
- Display: 3-digital LED
- Type of output: Constant voltage or Constant current
- Resume function: Yes, after power failure
- Operation conditions: 0-40°C, 0-95% humidity
- Safety features: No load detection, sudden load change detection, overload / short circuit protection.
- Dimension: 31×24×8 cm (L×W×H)
- Weight: 1.5kg
- Certificate: CE

Ordering Information

Cat. No.	Product Description
BL-PBC3220	High Current Power Supply 300V / 2000mA / 200W
BL-PBV6820	High Voltage Power Supply 600V / 800mA / 200W