




Goat Anti-Chicken IgY(H+L)

Cat. No.	Format	Size
6100-01	Purified (UNLB)	1.0 mg
6100-02	Fluorescein (FITC)	1.0 mg
6100-04	Alkaline Phosphatase (AP)	1.0 mL
6100-05	Horseradish Peroxidase (HRP)	1.0 mL
6100-08	Biotin (BIOT)	1.0 mg
6100-09	R-phycoerythrin (PE)	0.5 mg
6100-30	Alexa Fluor [®] 488 (AF488)	1.0 mg



ELISA plate was coated with purified chicken IgY and IgM. Immunoglobulins were detected with Goat Anti-Chicken IgY(H+L)-BIOT (SB Cat. No. 6100-08) followed by Streptavidin-HRP (SB Cat. No. 7100-05).

Description

Specificity	Reacts with the heavy and light chains of chicken IgY
Source	Pooled antisera from goats hyperimmunized with IgY
Cross Adsorption	None; may react with immunoglobulins from other species and the light chains of other chicken immunoglobulins
Purification	Affinity chromatography on chicken IgY covalently linked to agarose

Applications

Quality tested applications include –

ELISA ^{1,3,5-7}

FLISA

Other referenced applications include –

FC ⁸⁻¹⁰

IHC-PS ¹¹

ICC ^{4,5,12-16}

WB ^{1,3,4,17-23}

Working Dilutions

ELISA	AP conjugate	1:2,000 – 1:4,000
	HRP conjugate	1:4,000 – 1:8,000
	BIOT conjugate	1:5,000 – 1:20,000

FLISA	FITC and AF488 conjugates	1:200 – 1:400
	PE conjugate	≤ 1 µg/mL

Other Applications Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.

For Research Use Only. Not for Diagnostic or Therapeutic Use.

Handling and Storage

- The purified (UNLB) antibody is supplied as 1.0 mg purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 1.0 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- The alkaline phosphatase (AP) conjugate is supplied as 1.0 mL in a stock solution of 50 mM Tris/1 mM MgCl₂/50% glycerol, pH 8.0, containing NaN₃ as preservative. Store at 2-8°C or long-term at -20°C.
- The horseradish peroxidase (HRP) conjugate is supplied as 1.0 mL in a stock solution of 50% glycerol/50% PBS, pH 7.4. No preservative added. Store at 2-8°C or long-term at -20°C.
- The biotin (BIOT) conjugate is supplied as 1.0 mg in 2.0 mL of PBS/NaN₃. Store at 2-8°C.
- The R-phycoerythrin (PE) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- The Alexa Fluor[®] 488 (AF488) conjugate is supplied as 1.0 mg in 1.0 mL of PBS/NaN₃. Store at 2-8°C.
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

Warning

Some reagents contain sodium azide. Please refer to product specific SDS.

References

1. Chart H, Rowe B, Baskerville A, Humphrey TJ. Serological response of chickens to Salmonella enteritidis infection. *Epidemiol Infect.* 1990;104:63-71. (ELISA, WB)
2. Bouchard C, Galinha A, Tartour E, Fridman WH, Sautès C. A transforming growth factor β -like immunosuppressive factor in immunoglobulin G-binding factor. *J Exp Med.* 1995;182:1717-26. (ELISA)
3. Lundquist P, Ritchie HH, Moore K, Lundgren T, Linde A. Phosphate and calcium uptake by rat odontoblast-like MRPC-1 cells concomitant with mineralization. *J Bone Miner Res.* 2002;17:1801-13. (ELISA, WB)
4. Réfega S, Girard-Misguich F, Bourdieu C, Péry P, Labbé M. Gene discovery in Eimeria tenella by immunoscreening cDNA expression libraries of sporozoites and schizonts with chicken intestinal antibodies. *Vet Parasitol.* 2003;113:19-33. (WB, ICC)
5. Réfega S, Cluzeaud M, Péry P, Labbé M, Girard-Misguich F. Production of a functional chicken single-chain variable fragment antibody derived from caecal tonsils B lymphocytes against macrogamonts of Eimeria tenella. *Vet Immunol Immunopathol.* 2004;97:219-30. (ELISA, ICC)
6. Windisch M, Hess M. Establishing an indirect sandwich enzyme-linked-immunosorbent-assay (ELISA) for the detection of antibodies against Histomonas meleagridis from experimentally infected specific pathogen-free chickens and turkeys. *Vet Parasitol.* 2009;161:25-30. (ELISA)
7. van der Goot JA, Engel B, van de Water SG, Buist W, de Jong MC, Koch G, et al. Validation of diagnostic tests for detection of avian influenza in vaccinated chickens using Bayesian analysis. *Vaccine.* 2010;28:1771-7. (ELISA)
8. Lam KM. Chemotactic activities of avian lymphocytes. *Dev Comp Immunol.* 1999;23:641-7. (FC)
9. Sheela RR, Babu U, Mu J, Elankumaran S, Bautista DA, Raybourne RB, et al. Immune responses against Salmonella enterica serovar enteritidis infection in virally immunosuppressed chickens. *Clin Diagn Immunol.* 2003;10:670-9. (FC)
10. Nguyen DH, Tangvoranuntakul P, Varki A. Effects of natural human antibodies against a nonhuman sialic acid that metabolically incorporates into activated and malignant immune cells. *J Immunol.* 2005;175:228-36. (FC)
11. Spackman E, Day JM, Pantin-Jackwood MJ. Astrovirus, reovirus, and rotavirus concomitant infection causes decreased weight gain in broad-breasted white poults. *Avian Dis.* 2010;54:16-21. (IHC-PS)
12. Einfeld DA, Hunter E. Expression of the TM protein of Rous sarcoma virus in the absence of SU shows that this domain is capable of oligomerization and intracellular transport. *J Virol.* 1994;68:2513-20. (ICC)
13. Zekarias B, Mo H, Curtiss R 3rd. Recombinant attenuated Salmonella enterica serovar typhimurium expressing the carboxy-terminal domain of alpha toxin from Clostridium perfringens induces protective responses against necrotic enteritis in chickens. *Clin Vaccine Immunol.* 2008;15:805-16. (ICC)
14. Hu H, Roth JP, Estevez CN, Zsak L, Liu B, Yu Q. Generation and evaluation of a recombinant Newcastle disease virus expressing the glycoprotein (G) of avian metapneumovirus subgroup C as a bivalent vaccine in turkeys. *Vaccine.* 2011;29:8624-33. (ICC)
15. Li Y, Reddy K, Reid SM, Cox WJ, Brown IH, Britton P, et al. Recombinant herpesvirus of turkeys as a vector-based vaccine against highly pathogenic H7N1 avian influenza and Marek's disease. *Vaccine.* 2011;29:8257-66. (ICC)
16. Zhao W, Spatz S, Zhang Z, Wen G, Garcia M, Zsak L, et al. Newcastle disease virus (NDV) recombinants expressing infectious laryngotracheitis virus (ILTV) glycoproteins gB and gD protect chickens against ILTV and NDV challenges. *J Virol.* 2014;88:8397-406. (ICC)
17. Bernhisel-Broadbent J, Yolken RH, Sampson HA. Allergenicity of orally administered immunoglobulin preparations in food-allergic children. *Pediatrics.* 1991;87:208-14. (WB)
18. Larsen JK, Yamboliev IA, Weber LA, Gerthoffer WT. Phosphorylation of the 27-kDa heat shock protein via p38 MAP kinase and MAPKAP kinase in smooth muscle. *Am J Physiol.* 1997;273:L930-40. (WB)
19. McCracken AA, Werner ED, Powell MJ, Kruse KB, Brodsky JL. Differential fates of invertase mutants in the yeast endoplasmic reticulum. *Yeast.* 2000;16:49-55. (WB)
20. Zhang S, Xiao L, Zhou H, Yu Z, Chen H, Guo A, et al. Generation and characterization of an H5N1 avian influenza virus hemagglutinin glycoprotein pseudotyped lentivirus. *J Virol Methods.* 2008;154:99-103. (WB)
21. Küng E, Frey J. AvxA, a composite serine-protease-RTX toxin of Avibacterium paragallinarum. *Vet Microbiol.* 2013;163:290-8. (WB)
22. Zhang Z, Ma C, Zhao P, Duan L, Chen W, Zhang F, et al. Construction of recombinant marek's disease virus (rMDV) co-expressing AIV-H9N2-NA and NDV-F genes under control of MDV's own bi-directional promoter. *PLoS One.* 2014;9(3):e90677. (WB)
23. El-Ashram S, Sun X, Yin Q, Liu X, Suo X. Exploring early and late Toxoplasma gondii strain RH infection by two-dimensional immunoblots of chicken immunoglobulin G and M profiles. *PLoS One.* 2015;10(3):e0121647. (WB)

Alexa Fluor[®] 488, 647, and 555 are provided under an Intellectual property license from Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For information on purchasing a license to this product for any other use, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.