




## Mouse Anti-Hamster IgG<sub>1</sub>

Cat. No.	Format	Size
1940-01	Purified (UNLB)	0.5 mg
1940-02	Fluorescein (FITC)	0.5 mg
1940-05	Horseradish Peroxidase (HRP)	1.0 mL
1940-08	Biotin (BIOT)	0.5 mg
1940-09	R-phycoerythrin (PE)	0.1 mg



FLISA plate was coated with purified hamster IgG<sub>1</sub>, IgG<sub>2</sub>, and IgG<sub>3</sub>. Immunoglobulins were detected with serially diluted Mouse Anti-Hamster IgG<sub>1</sub>-FITC (SB Cat. No. 1940-02).

### Overview

<b>Clone</b>	SB139a
<b>Isotype</b>	Mouse (BALB/c) IgG <sub>2b</sub> K
<b>Immunogen</b>	Armenian hamster IgG
<b>Specificity</b>	Armenian/Syrian hamster IgG <sub>1</sub>

### Applications

ELISA – Quality tested <sup>1</sup>  
 FLISA – Quality tested  
 FC – Quality tested

### Working Dilutions

<b>ELISA</b>	Purified (UNLB) antibody	≤ 1 µg/mL
	HRP conjugate	1:4,000 – 1:8,000
	BIOT conjugate	1:5,000 – 1:10,000
<b>FLISA</b>	FITC conjugate	1:200 – 1:400
	PE conjugate	≤ 1 µg/mL
<b>Flow Cytometry</b>	FITC conjugate	≤ 1 µg/10 <sup>6</sup> cells
	PE conjugate	≤ 0.1 µg/10 <sup>6</sup> cells
	For flow cytometry, the suggested use of these reagents is in a final volume of 100 µL	
<b>Other Applications</b>	Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.	

**For Research Use Only. Not for Diagnostic or Therapeutic Use.**

## Handling and Storage

---

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of borate buffered saline, pH 8.2. *No preservatives or amine-containing buffer salts added.* Store at 2-8°C.
- The fluorescein (FITC) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN<sub>3</sub>. Store at 2-8°C.
- The horseradish peroxidase (HRP) conjugate is supplied as 1.0 mL of stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added. Store at 2-8°C or long-term at -20°C.
- The biotin (BIOT) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN<sub>3</sub>. Store at 2-8°C.
- The R-phycoerythrin (PE) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN<sub>3</sub> and a stabilizing agent. Store at 2-8°C. **Do not freeze!**
- Protect fluorochrome-conjugated forms from light. Reagents are stable for the period shown on the label if stored as directed.

## Warning

---

Some reagents contain sodium azide. Please refer to product specific (M)SDS.

## References

---

1. Ploquin A, Szécsi J, Mathieu C, Guillaume V, Barateau V, Ong KC, et al. Protection against henipavirus infection by use of recombinant adeno-associated virus-vector vaccines. J Infect Dis. 2013;207:469-78. (ELISA)