



Certificate of Analysis


Product:	Rabbit Monoclonal Antibody Biotin Anti-Human IgG2 Rabbit Monoclonal Antibody, Clone RM118
Catalog No.:	31-1020-02
Lot No.:	
Clone	RM118
Specificity	This antibody reacts to the heavy chain of human IgG2. No cross reactivity with human IgG1, IgG3, IgG4, IgM, IgA, IgD, IgE, mouse IgG, rat IgG, or goat IgG.
Application:	ELISA
Immunogen:	Human IgG2
Purity:	Protein A affinity purified from an animal origin-free culture supernatant
Size:	50 µg
Concentration:	1.0 mg/mL
Buffer:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Usage:	ELISA: 50ng/well – 200ng/well (for Capture); 0.05ug/mL – 0.2ug/mL (for Detection);
Storage and Stability:	Stable for 1 Year at -20.0°C from date of receipt.
Country of Origin:	U.S.A.
Intended Use:	For Research Use Only Not for Diagnostic or Therapeutic Use



Sandwich ELISA using RM118 as the capture antibody, and Biotinylated anti-human light chains ($\kappa+\lambda$) antibody RM129 as the detection antibody, followed by an AP conjugated streptavidin.



ELISA of human immunoglobulins shows RM118 only reacts to human IgG2. No cross reactivity with Human IgG1, IgG3, IgG4, IgE, IgD, IgA, mouse IgG, rat IgG, or goat IgG. The plate was coated with 50 ng/well of different immunoglobulins. A serial dilution of RM118 was used as the primary antibody. An alkaline phosphatase conjugated anti-rabbit IgG as the secondary antibody.



A titer ELISA of human IgG2. A serial dilution of RM118 was used as the primary antibody. An alkaline phosphatase conjugated anti-rabbit IgG as the secondary antibody.