



Certificate of Analysis

Product:	Rabbit Monoclonal Antibody Anti-Terminal deoxynucleotidyl transferase (TdT) Rabbit Monoclonal Antibody, Clone RM379
Catalog No.:	31-1265-00
Lot No.:	
Clone:	RM379
Specificity:	This antibody reacts to human Terminal deoxynucleotidyl transferase (TdT)
Application:	Immunohistochemistry, Western Blot
Immunogen:	A peptide corresponding to human Terminal deoxynucleotidyl transferase (TdT)
Purity:	Protein A affinity purified from an animal origin-free culture supernatant
Size:	100 µL
Buffer:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Usage:	IHC: 1:100 -1:200 dilution; WB: 1:1000 - 1:2000 dilution.
Storage and Stability:	Stable for 1 Year at -20.0°C from date of receipt.
Country of Origin:	U.S.A.
Intended Use:	For Research Use Only Not for Diagnostic or Therapeutic Use



Immunohistochemical staining of formalin fixed and paraffin embedded human thymus tissue section using anti-TdT rabbit monoclonal antibody (Clone RM379) at a 1:200 dilution.



Western Blot of Jurkat cells lysate using anti-TdT rabbit monoclonal antibody (Clone RM379) at a 1:2000 dilution.