



Certificate of Analysis

Product:	Rabbit Monoclonal Antibody Anti-Matrix metalloproteinase-12 (MMP-12) Rabbit Monoclonal Antibody, Clone RM381
Catalog No.:	31-1267-00
Lot No.:	
Clone:	RM381
Specificity:	This antibody reacts to human Matrix metalloproteinase-12 (MMP-12)
Application:	Immunohistochemistry, Western Blot
Immunogen:	A peptide corresponding to the C-terminus of human Matrix metalloproteinase-12 (MMP-12)
Purity:	Protein A affinity purified from an animal origin-free culture supernatant
Size:	100 μ L
Buffer:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Usage:	IHC: 1:1000 -1:2500 dilution; WB: 1:200 - 1:1000 dilution.
Storage and Stability:	Stable for 1 Year at -20.0°C from date of receipt.
Country of Origin:	U.S.A.
Intended Use:	For Research Use Only Not for Diagnostic or Therapeutic Use



Immunohistochemical staining of formalin fixed and paraffin embedded human lung tissue section using anti-MMP-12 rabbit monoclonal antibody (Clone RM381) at a 1:2500 dilution.



Western Blot of HEK 293 cells lysate using anti-MMP-12 rabbit monoclonal antibody (Clone RM381) at a 1:500 dilution.