


Certificate of Analysis

Product:	Rabbit Monoclonal Antibody Anti-DDDDK Tag Rabbit Monoclonal Antibody, Clone RM399
Catalog No.:	31-1285-00
Lot No.:	
Clone	RM399
Specificity	This antibody reacts to recombinant proteins containing DDDDK tag. No cross reactivity with other endogenous protein in mammalian or bacteria cells.
Application:	Western Blot, Immunoprecipitation,
Immunogen:	Mixture of a peptide with DDDDK-Tag at the N-terminal and another peptide with DDDDK-Tag at the C-terminal.
Purity:	Protein A affinity purified from an animal origin-free culture supernatant.
Size:	100 µg
Concentration:	1.0 mg/mL
Buffer:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Usage:	WB: 0.001 µg/mL – 0.1 µg/mL
Storage and Stability:	Stable for 1 Year at -20.0°C from date of receipt.
Country of Origin:	U.S.A.
Intended Use:	For Research Use Only Not for Diagnostic or Therapeutic Use



Western blot of 293T cells (1) transfected or (2) non-transfected with a DNA construct encoding DDDDK tag fusion protein, using anti-DDDDK tag rabbit monoclonal antibody, Clone RM399 at 0.001 µg/mL.