



Certificate of Analysis

Product:	Rabbit Monoclonal Antibody Anti-ZAP70 Rabbit Monoclonal Antibody, Clone RM408
Catalog No.:	31-1294-00
Lot No.:	
Clone:	RM408
Specificity:	This antibody reacts to human ZAP70. This antibody may also react to mouse or rat ZAP70, as predicted by immunogen homology.
Application:	Immunohistochemistry, Western Blot
Immunogen:	A peptide corresponding to the N-terminus of human ZAP70
Purity:	Protein A affinity purified from an animal origin-free culture supernatant
Size:	100 µL
Buffer:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Usage:	IHC: 1:100 -1:200 dilution WB: 1:500 - 1:1000 dilution.
Storage and Stability:	Stable for 1 Year at -20.0°C from date of receipt.
Country of Origin:	U.S.A.
Intended Use:	For Research Use Only Not for Diagnostic or Therapeutic Use



Immunohistochemical staining of formalin fixed and paraffin embedded human tonsil tissue section using Anti-ZAP70 Rabbit Monoclonal Antibody (Clone RM408) at a 1:200 dilution.



Western Blot of Jurkat cell lysate using Anti-ZAP70 rabbit monoclonal antibody (Clone RM408) at a 1:500 dilution.