


## Vertical Electrophoresis (BL- CZ-B)



### Technical specifications:

Glass plate size (W×L)	216×110 mm
Gel Size (W×L)	186×95 mm
Gel thickness	1.0(mm)
Gel amount	1~2 (pieces)
Sample volume	(1.0mm) 25、40、52wells
Buffer Volume	1800 (ml)
Dimension (L×W×H)	280×140×135 (mm)
Net weight	3.2 (kg)

### Characteristic features:

- The installation of the glass plate in just 15 seconds, and does not need the knob, fast and convenient.
- With the gel making device of normal position don't should to moving the glass plate from the process of making gel to electrophoresis.
- With the Precast gel inner core, do not affecting the operation of gel making while the electrophoresis.
- Safety cover open button design, convenient opening the cover.
- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Auto-switch-off when the lid is removed.