


Vertical Electrophoresis (BL- JX5L)



Technical specifications:

Glass plate size (W×L)	216×110 mm
Gel Size (W×L)	186×205 mm
Gel thickness	1.0(mm)
Gel amount	1~2 (pieces)
Sample volume	(1.0mm) 25、40、52wells
Buffer Volume	3500 (ml)
Dimension (L×W×H)	280×140×250 (mm)
Net weight	5.6 (kg)

Characteristic features:

- The installation of the glass plate in just 15 seconds, and does not need the knob, fast and convenient.
- With the gel making device modular, you can make multi slab gel at the same time.
- Safety cover open button design, convenient opening the cover.
- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Auto-switch-off when the lid is removed.