

Vertical Electrophoresis (BL- SCZ2+)



Technical specifications:

Glass plate size (W×L)	100×100 mm
Gel Size (W×L)	82×88 mm
Gel thickness	0.75、1.0、1.5 (mm)
Gel amount	1~2 (pieces)
Sample volume	(0.75mm) 11、15wells: (1.0mm) 11、15wells ; (1.5mm) 11、15wells
Buffer Volume	750 (ml)
Dimension (L×W×H)	150×120×115 (mm)
Net weight	1.5 (kg)

Characteristic features:

- The installation of the glass plate in just 15 seconds, and does not need the knob, fast and convenient.
- With a gel casting in original position, no need to move glasses from gel making to electrophoresis running. It is more convenient to observe gel preparing through both sides of the glasses.
- High transparent polycarbonate material injection molding, impact resistance, high temperature resistance, corrosion resistance.
- Safety cover open button design, convenient opening the cover.
- With the background color, easy to be observed in the sample adding and electrophoresis process.
- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Auto-switch-off when the lid is removed.